



Gabriel, CeraPlus™ Product User

Your Skin Health & Security Accessories Guide

Customise your ostomy pouching
system with confidence



CeraRing™ Seals

The confidence of a secure seal to help prevent leakage.

- Infused with ceramide to protect your skin
- Fill dips and creases around your stoma, for a smooth surface to improve barrier fit
- Stretch, mould or stack to customise with confidence
- Alcohol free, no sting



CeraRing™ Slim & Standard Seals

Hollister Code	Size	SAS Issue
Slim 8815	<i>Thickness</i> 2.3mm	3 boxes per month (10 per box)
Standard 8805	4.5mm	



CeraRing™ Round Convex Seals

Hollister Code	Size	SAS Issue
89520	<i>Inner Diameter</i> 20-25mm	3 boxes per month (10 per box)
89530	30-35mm	
89540	40-45mm	



CeraRing™ Oval Convex Seals

Hollister Code	Size	SAS Issue
89601	22×38mm	3 boxes per month (10 per box)
89602	30×48mm	
89603	38×56mm	



Scan for more information or to order a sample of CeraRing™ Slim or Standard Seal

Using CeraRing™ Seals

Option 1:

Apply flat or convex CeraRing™ seal to the adhesive side of your skin barrier.



Stretch and shape seal for a customised fit.

Option 2:

Apply flat or convex CeraRing™ seal to your skin, then apply your pouching system. Make sure the seal fits where the skin and stoma meet.




Scan for more information or to order a sample of CeraRing™ Round Convex or Oval Convex Seal

Adapt Barrier Extenders

The confidence of enhancing your pouch security.

Hollister Code	SAS Issue	Product Info
79402	2 boxes per month (30 per box)	Customise usage with 3 separate sections Hydrocolloid material



Enhance the feeling of pouch security by adding a frame to your skin barrier, also helping to extend pouch wear time.



Increase the feeling of security when undertaking **physical activity**



Suitable for use in **warm climates** to reduce barrier edge lifting

Using Adapt Barrier Extenders

1. Remove existing skin barrier and clean the skin around the stoma. Apply new skin barrier.
2. Remove paper backing by stretching lightly at each end of the extender.
3. Release only one side of the paper backing.
4. Place the adhesive end of the barrier extender along the perimeter of the skin barrier; covering both skin barrier and the skin next to it. Remove remaining paper backing.
5. Repeat for any additional strips you require around your barrier.
6. Barrier extenders will remove easily with your routine pouch change.



Scan for more information or to order a sample of Barrier Extenders

Adapt Remover Wipes and Spray

The confidence of an easy pouch removal that is gentle on skin.

Adapt No-Sting Adhesive Remover Wipes

Hollister Code	SAS Issue	Product Info
7760	14 boxes per year (50 sachets per box)	Alcohol free. No sting

Adapt No-Sting Adhesive Remover Spray

Hollister Code	SAS Issue	Product Info
7737	2 spray cans per month	Alcohol free. No sting

Using Adhesive Remover Spray

1. Shake well before use.
2. Peel back one edge of skin barrier.
3. Hold can so arrow faces the base of the skin barrier, and spray until the adhesive releases.
4. Ensure skin is dry and apply new pouching system as usual.

Spray

Spray until the skin barrier releases.



- Gently remove adhesive products and residue to make removal of skin barriers fast and easy
- Help reduce potential for skin stripping



Using Adhesive Remover Wipes

1. Press or rub wipe against your skin as you peel off your skin barrier and/or tape.
2. Press or rub wipe over any areas of the skin where tape residue remains.
3. Ensure skin is dry and apply new ostomy pouching system as usual.

Wipes

Rub the wipe over skin as you peel off skin barrier.



Wipes

Clean off remaining tape residue.



Scan for more information or to order a sample of Adhesive Remover Wipes

Adapt Lubricating Deodorant

The confidence of eliminating odour and easy pouch emptying.

Hollister Code	SAS Issue	Product Info
78500 (Bottle)	1 bottle per month	All in one solution; lubricant and deodorant
78501 (Sachets)	2 sachet boxes per year (50 sachets per box)	in same bottle

Using for Drainable Pouches

1. Apply pouch.
2. Bottle: add approximately 5ml or 1tsp (Figure 1).
Sachet: add the complete contents of the packet (Figure 2).
Close the pouch in the normal manner.
3. Rub to thoroughly coat the inside of the pouch (Figure 3).
Repeat Steps 1 or 2 after every emptying.

Using for Closed Pouches

1. Do not remove release paper on barrier to add product.
2. Bottle: add approximately 5ml or 1tsp (Figure 4).
Sachet: add the complete contents of the packet.
3. Apply your pouch. Rub to thoroughly coat the inside of the pouch (Figure 3).



Lubricate inside of the pouch to limit pancaking and for easy emptying.



Scan for more information or to order a sample of Lubricating Deodorant

“My nurse recommended using CeraPlus soft convex skin barrier with the CeraRing, which are both infused with ceramide ... my skin feels healthy and looks amazing!”

SHAUN, CERAPLUS™ PRODUCT USER



“I live a very active lifestyle, and barrier extenders give me the security and confidence to get on with my life without worry.”

LATESHA, CERAPLUS™ PRODUCT USER



The testimonials, statements, and opinions presented are applicable to the people depicted. These testimonials are representative of their experiences, but the exact results and experience will be unique and individual to each person. LaTasha and Shaun received compensation from Hollister Inc. for this testimonial.

Skin Health and Security Ostomy Accessories

Hollister Code	Size	SAS Issue <i>(Amount of product you are able to order through the Stoma Appliance Scheme)</i>
----------------	------	--



CeraRing Slim & Standard Seals

	<i>Thickness</i>	
Slim 8815	2.3mm	3 boxes per month (10 per box)
Standard 8805	4.5mm	



CeraRing Round Convex Seals

	<i>Inner Diameter</i>	
89520	20-25mm	3 boxes per month (10 per box)
89530	30-35mm	
89540	40-45mm	



CeraRing Oval Convex Seals

89601	22×38mm	3 boxes per month (10 per box)
89602	30×48mm	
89603	38×56mm	



Adapt Barrier Extenders

79402	3 sections	2 boxes per month (30 per box)
-------	------------	--------------------------------



Adapt No-Sting Adhesive Remover Wipes

7760	50 sachets per box	14 boxes per year
------	--------------------	-------------------

Adapt No-Sting Adhesive Remover Spray

7737	50ml	2 spray cans per month
------	------	------------------------



Adapt Lubricating Deodorant

78500 (Bottle)	236mL bottle	1 bottle per month
78501 (Sachets)	50 sachets per box	2 sachet boxes per year

Ready to personalise your pouching system? We're here to help.
Customer Care: 1800 880 851

Your Stomal Therapy Nurse is able to help guide you if you are unsure on which accessories are right for you.

Prior to use, be sure to read the instructions for use for information regarding Intended Use, Contraindications, Warnings, Precautions, and Instructions for Use.