



Pouchkins

Paediatric Pouching Systems

**With the gentleness and flexibility of SoftFlex**

- Designed to be gentle on sensitive infant skin
- No starter hole which enables the opening to be placed where it is needed
- No tape by design
- Flexible and sized for babies or children



# Pouchkins

## Paediatric Pouching Systems



	Newborn	Infant	Infant	Child
	3778	3796	3797	3795
	Two Piece Drainable	One Piece Drainable	One Piece Urostomy	One Piece Drainable
<b>Skin Barrier</b> Cut-To-Fit (CTF) Flat, without Tape Border, Opening Up to;	35 x 22 mm	38 mm	38 mm	51 mm
<b>Gentle &amp; Flexible SoftFlex Skin Barrier</b> (Without starter hole for versatility)	✓	✓	✓	✓
<b>Pouch Length</b>	16.5 cm	16 cm	22 cm	18 cm
<b>Pouch Volume</b>	178 ml	193 ml	225 ml	274 ml
<b>Odour-Barrier Pouch Film</b>	✓	✓	✓	✓
<b>Colour</b>	Ultra Clear	Transparent	Transparent	Ultra Clear
<b>Soft Backing</b>		✓	✓	✓
<b>Filter</b>				✓
<b>Lock n' Roll Closure Microseal Closure</b>		✓		✓
<b>Box QTY</b>	15	10	10	10
	<ul style="list-style-type: none"> <li>Also ideal for fistula management.</li> <li>Drain tap at bottom, compatible with luer lock syringe.</li> <li>Secure drainable cap.</li> <li>Can be used with or without the drain valve, cut at the bottom to use as drainable pouch with soft wire closure</li> </ul>	<ul style="list-style-type: none"> <li>Non-woven ComfortWear panels on the body side for enhanced comfort and helps prevent the pouch from sticking to the skin.</li> </ul>	<ul style="list-style-type: none"> <li>Anti-reflux valve.</li> <li>Wide bore tubing for easy drainage.</li> <li>Secure drainable cap.</li> </ul>	<ul style="list-style-type: none"> <li>With integrated AF300 Filter allowing for slow release and deodorisation of gas.</li> </ul>

- ✓ All Pouchkins paediatric products use **SoftFlex** - a thin, gentle and flexible standard wear skin barrier. The formulation of **SoftFlex** skin barrier incorporates cotton fibres to enhance strength and integrity.



**Hollister Australia**  
PO Box 375  
Box Hill, VIC 3128  
1800 880 840  
[www.hollister.com.au](http://www.hollister.com.au)

**Hollister New Zealand**  
PO Box 62-027  
Mt Wellington, Auckland  
0800 678 669  
[www.hollister.co.nz](http://www.hollister.co.nz)