

**Hollister Continence Care**

Apogee Intermittent Catheter

---

People First.



Trusted Quality.  
Reliable Performance.  
Cost Effective.


## Apogee Intermittent Catheter

Flexibility and control - manage your continence needs with confidence.  
Designed for comfort with features you need:

- **Non-lubricated**
- **Affordable** – meeting the needs of all ISC users
- **Smooth eyelets** provide comfort on insertion and removal
- **Smooth** 100% medical grade PVC - Latex free
- **Easy to open** packaging
- **Full range of sizes** from 10-18Fr
- **Male and Female lengths**
- **Firm but flexible** for greater comfort and control

## The eyelets on all Hollister catheters are smooth and rounded.

Compare ours with other brands and you will feel the difference. Formed with ultrasonic technology, these catheter eyelets provide maximum comfort on insertion and removal, and help to minimise trauma to the urethra.





# Ordering Information

<b>Apogee: Intermittent Catheter</b>						
<b>Hollister Stock Code</b>	<b>Colour</b>	<b>Description</b>	<b>Gender</b>	<b>Length</b>	<b>Size</b>	<b>Box QTY</b>
11006	●	Apogee IC Non-lubricated	Female	15 cm	10Fr	30
11206	○	Apogee IC Non-lubricated	Female	15 cm	12Fr	30
11406	●	Apogee IC Non-lubricated	Female	15 cm	14Fr	30
11016	●	Apogee IC Non-lubricated	Male	40 cm	10Fr	30
11216	○	Apogee IC Non-lubricated	Male	40 cm	12Fr	30
11416	●	Apogee IC Non-lubricated	Male	40 cm	14Fr	30
11616	●	Apogee IC Non-lubricated	Male	40 cm	16Fr	30
11816	●	Apogee IC Non-lubricated	Male	40 cm	18Fr	30
10616	●	Apogee IC Non-lubricated	Male	40 cm	6 Fr	30
10816	●	Apogee IC Non-lubricated	Male	40 cm	8 Fr	30
11226	○	Apogee IC Non-lubricated Tiemann Tip	Male	40 cm	12 Fr	30
11426	●	Apogee IC Non-lubricated Tiemann Tip	Male	40 cm	14 Fr	30

**For more information call  
Customer Service on  
0800 678 669**

**[www.hollister.co.nz](http://www.hollister.co.nz)**

