



Freedom to Choose
Freedom to Live

2020/21 Product & Resource Catalogue

Territory Manager

Contacts

DANSAC

Gene Wotherspoon

Mobile: 021 886 133

Email: gene.wotherspoon@dansac.co.nz

HOLLISTER

Mark Liddle

Mobile: 021 907 587

Email: mark.liddle@hollister.co.nz

CUSTOMER SERVICE

Freecall: 0800 678 669

Email: NZcustomerservice@libmed.co.nz
gale.wenski@libmed.co.nz

LIBERTY MEDICAL NEW ZEALAND

58 Richard Pearse Drive
Airport Oaks
Auckland
New Zealand

Contents

One-Piece Pouches	5	Paediatric Pouches	55
Closed Flat	8	Accessories	59
Closed Soft Convex	12	Seals	62
Closed Firm Convex	14	Skin Care	64
Drainable Flat	16	Irrigation Products	68
Drainable Soft Convex	20	Intermittent Catheter	
Drainable Firm Convex	24	& Drainage Bags	71
Urostomy Flat	28	Tube Attachment	
Urostomy Soft Convex	30	Devices (TAD)	75
Urostomy Firm Convex	32	Resources	81
Two-Piece Barriers	35	Index	101
Flat	38		
Convex	42		
Two-Piece Pouches	45		
Closed	48		
Drainable	50		
Urostomy	52		

ONE-PIECE POUCHES

Liberty
freedom to live

The CeraPlus™ Skin Barrier with Remois Technology*

The importance of peristomal skin health in ostomy care can't be overstated. The condition of this skin greatly affects how well the skin barrier can adhere, as well as the general comfort and well-being of a person living with a stoma.

The Ceramide Advantage

Ceramide is a natural component of human skin, consisting of waxy lipid molecules. Ceramide links cells of the outer skin (epidermis) together to form a waterproof, protective barrier. Ceramides help prevent water loss that can lead to skin damage and dryness.

The CeraPlus Skin Barrier is Designed to Support Healthy Skin From the Start

It features a formulation infused with ceramide to help:

- Protect the skin's natural moisture barrier
- Helps keep your peristomal skin healthy
- Helps minimise moisture loss from the skin (transepidermal water loss – TEWL).



The Science of Skin Health Meets the Ultimate Urostomy Pouch.




* Remois is a technology of Alcare Co., Ltd.



Dansac TRE™

Technology Explained

- The best skin is healthy skin.
- pH is a measure of how acidic or alkaline something is. pH 7 is neutral.
- Healthy skin has a pH between 4 and 6 which means it is slightly acidic. This is known as the acid mantle.
- Maintaining the acid mantle of the skin is essential to help protect its antimicrobial properties and maintain skin integrity.
- The acid mantle of the skin may become weakened if stoma output (which includes digestive enzymes), or perspiration stay in contact with the skin.
- Dansac's TRE technology is designed to protect and maintain the skin's acid mantle, so that it can go on working to protect you.



Moderma Flex CeraPlus Pouch

Box of 30

- **Ceramide** Infused Barrier
- AF300 Filter



**VIEWING
OPTION**

MIDI		
CLEAR	VIEWING OPTION	STARTER HOLE
	52100	15 – 55 mm (CTF)
	52125	25 mm
	52130	30 mm
	52135	35 mm

Moderma Flex Pouch

Box of 30

- **SoftFlex** Skin Barrier
- AF300 Filter
- Dimples Technology



MIDI		
CLEAR	OPAQUE	STARTER HOLE
22500	22100	15 – 55 mm (CTF)
	22120	20 mm
	22125	25 mm
	22130	30 mm
	22135	35 mm
	22140	40 mm

MAXI		
CLEAR	OPAQUE	STARTER HOLE
22400	22300	15 – 55 mm (CTF)
	22325	25 mm
	22330	30 mm
	22335	35 mm

Moderma Flex QuietWear Pouch

Box of 30


- **SoftFlex** Skin Barrier
- AF300 Filter
- Dimples Technology



MIDI		
CLEAR	OPAQUE	STARTER HOLE
	42100	15 – 55 mm (CTF)
	42125	25 mm
	42130	30 mm
	42135	35 mm
	42140	40 mm
	42145	45 mm

NovaLife TRE 1 Closed Flat Pouch

Box of 30



MIDI		
CLEAR	EASIVIEW	STARTER HOLE
	3001-15	15 – 55/70 mm (CTF)
	3001-25	25 – 55/70 mm (PC)
	3001-30	30 – 55/70 mm (PC)
	3001-35	35 – 55/70 mm (PC)
	3001-40	40 – 55/70 mm (PC)
	3001-45	45 – 55/70 mm (PC)

MAXI		
CLEAR	EASIVIEW	STARTER HOLE
	3007-15*	15 – 55/70 mm (CTF)
	3005-15	15 – 70/90 mm (CTF)
	3007-25*	25 – 55/70 mm (PC)
	3007-30*	30 – 55/70 mm (PC)
	3005-25	25 – 70/90 mm (PC)
	3005-30	30 – 70/90 mm (PC)

* MIDI barrier on MAXI pouch

NovaLife GX 1 Closed Flat Pouch

Box of 30



MINI		
CLEAR	OPAQUE	STARTER HOLE
	903-30	30 – 38/50 mm (PC)
	903-35	35 – 38/50 mm (PC)

MIDI		
CLEAR	EASIVIEW	STARTER HOLE
902-10	9010-10	10 – 55/70 mm (CTF)
902-25	9010-25	25 – 55/70 mm (PC)
902-30	9010-30	30 – 55/70 mm (PC)
902-35	9010-35	35 – 55/70 mm (PC)
902-40	9010-40	40 – 55/70 mm (PC)
902-45	9010-45	45 – 55/70 mm (PC)
902-50	9010-50	50 – 55/70 mm (PC)

MAXI		
CLEAR	EASIVIEW	STARTER HOLE
906-10	9050-10	10 – 70/90 mm (CTF)
	9050-25	25 – 70/90 mm (PC)
	9050-30	30 – 70/90 mm (PC)
	9050-35	35 – 70/90 mm (PC)

All pre-cut pouches can still be cut out to listed size (ie- 55/70 mm)

Moderna Flex Oval Pouch

Box of 30

- **SoftFlex** Skin Barrier
- AF300 Filter

MAXI		
CLEAR	OPAQUE	STARTER HOLE
22401	22301	15 – 75 mm (CTF)



Stoma Cap

Box of 30

- **SoftFlex** Skin Barrier

MINI		
CLEAR	OPAQUE	STARTER HOLE
	3184	51 mm



Compact Stoma Cap

Box of 30

- **SoftFlex** Skin Barrier

MINI		
CLEAR	OPAQUE	STARTER HOLE
	3192	30 mm
	3193	40 mm



Nova 1 Closed Flat Pouch

Box of 30



MINI/STOMA CAP		
CLEAR	OPAQUE	STARTER HOLE
	829-20	20 – 50 mm (CTF)
	829-25	25 mm (PC)
	829-30	30 mm (PC)
	829-40	40 mm (PC)

MINI		
CLEAR	OPAQUE	STARTER HOLE
	803-20	20 – 50 mm (CTF)
	803-25	25 mm (PC)
	803-30	30 mm (PC)
	803-35	35 mm (PC)

MIDI		
CLEAR	OPAQUE	STARTER HOLE
802-20	801-20	20 – 60 mm (CTF)
802-25	801-25	25 mm (PC)
802-30	801-30	30 mm (PC)
802-35	801-35	35 mm (PC)
	801-40	40 mm (PC)
	801-45	45 mm (PC)
	801-50	50 mm (PC)

MAXI		
CLEAR	OPAQUE	STARTER HOLE
806-20	805-20	20 – 60 mm (CTF)

Moderma Flex Pouch

- Box of 30
- **SoftFlex** Skin Barrier
 - Soft Convexity
 - Belt Tabs
 - AF300 Filter
 - Convexity Depth = 6 mm



VIEWING
OPTION

MIDI		
CLEAR	VIEWING OPTION	STARTER HOLE
	2261038	15 – 38 mm (CTF)
	2261055	15 – 55 mm (CTF)
	2261125	25 mm
	2261130	30 mm
	2261135	35 mm
	2261140	40 mm

MAXI		
CLEAR	VIEWING OPTION	STARTER HOLE
	2281025	15 – 25 mm (CTF)
	2281038	15 – 38 mm (CTF)
	2281055	15 – 55 mm (CTF)
	2281125	25 mm
	2281130	30 mm
	2281135	35 mm
	2281140	40 mm

NovaLife TRE 1 Closed Soft Convex Pouch

Box of 10



MIDI		
CLEAR	EASIVIEW	STARTER HOLE
	3071-24	15 – 24 mm (CTF)
	3071-34	15 – 34 mm (CTF)
	3071-25	25 mm (PC)
	3071-30	30 mm (PC)
	3071-35	35 mm (PC)

MAXI		
CLEAR	EASIVIEW	STARTER HOLE
	3071-44	15 – 44 mm (CTF)
	3071-54	15 – 54 mm (CTF)
	3071-40	40 mm (PC)

NovaLife GX 1 Closed Soft Convex Pouch

Box of 10



MIDI		
CLEAR	EASIVIEW	STARTER HOLE
	971-34	15 – 35 mm (CTF)
	971-25	25 mm (PC)
	971-30	30 mm (PC)
	971-35	35 mm (PC)

MAXI		
CLEAR	EASIVIEW	STARTER HOLE
	971-44	15 – 44 mm (CTF)
	971-54	15 – 55 mm (CTF)
	971-40	40 mm (PC)

Nova 1 Closed Soft Mouldable Pouch

Box of 10



MIDI		
CLEAR	OPAQUE	STARTER HOLE
	871-25	25 – 35 mm (stretch)
	871-35	35 – 45 mm (stretch)

Moderma Flex CeraPlus Pouch

Box of 30

- **Ceramide** Infused Barrier
- AF300 Filter
- Belt Tabs
- Convexity Depth = 4.16 mm



**VIEWING
OPTION**

MIDI		
CLEAR	VIEWING OPTION	STARTER HOLE
	52600	15 – 38 mm (CTF)
	52601	15 – 51 mm (CTF)
	52625	25 mm
	52630	30 mm
	52635	35 mm

Moderma Flex Pouch

Box of 30

- **SoftFlex** Skin Barrier
- Cushion Fit
- Belt Tabs
- AF300 Filter
- Convexity Depth = 4.16 mm



MAXI		
CLEAR	OPAQUE	STARTER HOLE
22900	22800	15 – 38 mm (CTF)
22901	22801	15 – 51 mm (CTF)
	22820	20 mm
22925	22825	25 mm
22930	22830	30 mm
22935	22835	35 mm

NovaLife TRE 1 Closed Convex Pouch

Box of 10



NEW

MIDI		
CLEAR	EASIVIEW	STARTER HOLE
	3031-37	15 – 37 mm (CTF)
	3031-59	15 – 59 mm (CTF)
	3035-25	25 mm (PC)
	3031-30	30 mm (PC)
	3031-35	35 mm (PC)

NovaLife GX 1 Closed Convex Pouch

Box of 10



MIDI		
CLEAR	EASIVIEW	STARTER HOLE
	9310-37	10 – 37 mm (CTF)
932-25	9310-25	25 mm (PC)
932-30	9310-30	30 mm (PC)
932-35	9310-35	35 mm (PC)
932-40		40 mm (PC)
932-45	9310-45	45 mm (PC)

MAXI		
CLEAR	EASIVIEW	STARTER HOLE
	9310-50	50 – 59 mm (CTF)

Nova 1 Closed Convex Pouch

Box of 10



MIDI		
CLEAR	OPAQUE	STARTER HOLE
	831-25	25 mm (PC)

MAXI		
CLEAR	OPAQUE	STARTER HOLE
832-37	831-37	15 – 37 mm (CTF)
	831-46	15 – 46 mm (CTF)
	831-30	30 mm (PC)
	831-35	35 mm (PC)
	831-40	40 mm (PC)
	831-45	45 mm (PC)

Moderma Flex CeraPlus Pouch – with Tape Border

Box of 30

- **Ceramide** Infused Barrier
- Lock 'n Roll Microseal Closure
- AF300 Filter



**VIEWING
OPTION**

MIDI		
CLEAR	VIEWING OPTION	STARTER HOLE
53500	54200	15 – 64 mm (CTF)
MAXI		
CLEAR	VIEWING OPTION	STARTER HOLE
53400	54300	15 – 64 mm (CTF)
	54325	25 mm
	54330	30 mm
	54335	35 mm

Moderma Flex FlexWear Pouch – with Tape Border

Box of 30

- **FlexWear** Durable Skin Barrier
- Lock 'n Roll Microseal Closure
- AF300 Filter



MAXI		
CLEAR	OPAQUE	STARTER HOLE
58400	58300	15 – 64 mm (CTF)
	58325	25 mm

Moderma Flex Pouch

Box of 10

- **SoftFlex** Skin Barrier
- Lock 'n Roll Microseal Closure
- AF300 Filter



MAXI		
CLEAR	OPAQUE	STARTER HOLE
28400	28300	15 – 55 mm (CTF)
	28325	25 mm
	28330	30 mm
	28335	35 mm
	28340	40 mm

NovaLife TRE 1 Open Flat Pouch

Box of 30 (Mini/Midi)

Box of 10 (Maxi)



NEW

MINI		
CLEAR	OPAQUE	STARTER HOLE
3026-15	3025-15	15 – 38/50 mm (CTF)
MIDI		
CLEAR	EASIVIEW	STARTER HOLE
3024-15	3023-15	15 – 55/70 mm (CTF)
	3023-25	25 – 55/70 mm (PC)
	3023-30	30 – 55/70 mm (PC)
	3023-35	35 – 55/70 mm (PC)
	3023-40	40 – 55/70 mm (PC)
MAXI		
CLEAR	EASIVIEW	STARTER HOLE
	3017-15*	15 – 55/70 mm (CTF)
3016-15	3015-15	15 – 70/90 mm (CTF)
	3017-25*	25 – 55/70 mm (PC)
	3017-30*	30 – 55/70 mm (PC)
	3015-25	25 – 70/90 mm (PC)
	3015-30	30 – 70/90 mm (PC)

* MIDI barrier on MAXI pouch

NovaLife GX 1 Open Flat Pouch

Box of 30 (Mini/Midi)

Box of 10 (Maxi)



MINI		
CLEAR	OPAQUE	STARTER HOLE
	925-10	10 – 38/50 mm (CTF)
MIDI		
CLEAR	OPAQUE	STARTER HOLE
924-10	923-10	10 – 55/70 mm (CTF)
	923-25	25 – 55/70 mm (PC)
	923-30	30 – 55/70 mm (PC)
	923-35	35 – 55/70 mm (PC)
MAXI		
CLEAR	OPAQUE	STARTER HOLE
916-10	915-10	25 – 70/90 mm (CTF)

All pre-cut pouches can still be cut out to listed size (ie- 55/70 mm)

[Contents](#)

[Index](#)

Moderma Flex Oval Pouch

Box of 30

- **SoftFlex** Skin Barrier
- Lock 'n Roll Microseal Closure
- AF300 Filter

MAXI		
CLEAR	OPAQUE	STARTER HOLE
28401	28301	15 – 75 mm (CTF)



Post-Operative Pouch

Box of 5

- Premier Post-op
- Lock 'n Roll Microseal Closure
- No Filter
- Ultra Clear Film

MAXI		
CLEAR	OPAQUE	STARTER HOLE
89004		12 – 64 mm (CTF)



Moderma Flex High Output Pouch

Box of 10

- **Flextend** Durable Skin Barrier
- Soft Tap
- Tape Border
- AF300 Filter

MAXI		
CLEAR	OPAQUE	STARTER HOLE
26064		15 – 64 mm (CTF)



NovaLife GX+ 1 Open Flat Pouch

Box of 30 (Midi)

Box of 10 (Maxi)



MIDI		
CLEAR	OPAQUE	STARTER HOLE
	9630-25	25 – 55/70 mm (PC)
	9630-30	30 – 55/70 mm (PC)
	9630-35	35 – 55/70 mm (PC)

MAXI		
CLEAR	OPAQUE	STARTER HOLE
	9610-10	10 – 70/90 mm (CTF)
	9610-25	25 – 70/90 mm (PC)

Nova 1 FoldUp Flat Pouch

Box of 30 (Mini/Midi)

Box of 10 (Maxi)



MINI		
CLEAR	OPAQUE	STARTER HOLE
826-15*	825-15	15 – 50 mm (CTF)

MIDI		
CLEAR	OPAQUE	STARTER HOLE
824-15*	823-15	15 – 60 mm (CTF)
	823-25	25 mm (PC)
	823-30	30 mm (PC)
	823-35	35 mm (PC)
	823-40	40 mm (PC)

MAXI		
CLEAR	OPAQUE	STARTER HOLE
816-15*	815-15	15 – 90 mm (CTF)

Nova 1 FoldUp X3 Pouch

Box of 10



MIDI		
CLEAR	OPAQUE	STARTER HOLE
	861-15	15 – 46 mm (CTF)
	861-30	30 mm (PC)

Moderma Flex CeraPlus Pouch – with Tape Border

- Box of 30
- **Ceramide** Infused Barrier
 - Lock 'n Roll Closure
 - AF300 Filter
 - Convexity Depth = 6 mm



VIEWING
OPTION

MAXI		
CLEAR	VIEWING OPTION	STARTER HOLE
53100	54100	15 – 38 mm (CTF)
53101	54101	15 – 55 mm (CTF)
	54125	25 mm
	54130	30 mm
	54135	35 mm

Moderma Flex Pouch

- Box of 10
- **FlexWear** Skin Barrier
 - Lock 'n Roll Closure
 - AF300 Filter
 - Convexity Depth = 6 mm



MIDI		
CLEAR	VIEWING OPTION	STARTER HOLE
	5681038	15 – 38 mm (CTF)
	5681055	15 – 51 mm (CTF)
	5681120	20 mm
	5681125	25 mm
	5681130	30 mm
	5681135	35 mm

NovaLife TRE 1 Open Soft Convex Pouch

Box of 10



MIDI		
CLEAR	EASIVIEW	STARTER HOLE
	3081-24	15 – 24 mm (CTF)
	3081-34	15 – 34 mm (CTF)
	3081-25	25 mm (PC)
	3081-30	30 mm (PC)
	3081-35	35 mm (PC)

MAXI		
CLEAR	EASIVIEW	STARTER HOLE
3084-34	3083-34	15 – 34 mm (CTF)
	3081-44	15 – 44 mm (CTF)
	3081-54	15 – 54 mm (CTF)
	3083-25	25 mm (PC)
	3083-30	30 mm (PC)
	3083-35	35 mm (PC)
	3081-40	40 mm (PC)
	3081-45	45 mm (PC)

NovaLife GX 1 Open Soft Convex Pouch

Box of 10



MIDI	
EASIVIEW	STARTER HOLE
981-24	15 – 24 mm (CTF)
981-34	15 – 34 mm (CTF)
981-25	25 mm (PC)
981-30	30 mm (PC)
981-35	35 mm (PC)

MAXI	
EASIVIEW	STARTER HOLE
983-24	15 – 24 mm (CTF)
983-34	15 – 34 mm (CTF)
981-44	15 – 44 mm (CTF)
981-54	15 – 54 mm (CTF)
981-40	40 mm (PC)
981-45	45 mm (PC)
983-30	30 mm (PC)
983-35	35 mm (PC)



AF300 Filter – A Filter You Can Use With Confidence

- May help reduce gas build up and ballooning.
- Internal and external liquid protection. No stickers/ covers required.



Tape Border – Designed With Security In Mind

- Breathable
- Non-woven adhesive (not a microporous adhesive)
- Lightweight and highly flexible
- Conformable
- Creates a lower barrier profile



Viewing Option

The discretion of beige, and the confidence of clear – all in one pouch.


- The viewing option is utilised by breaking the perforated seal.
- When the flap is lifted, the stoma and output can be observed clearly and privately.
- The curved seal helps prevent water from collecting during bathing or showering.



Security Made Simple

Newly designed security flap and tail reinforcement.

- Lock 'n Roll Microseal Closure
- Transparent security flap – easy to open and close
- Flexible ribbed strips reinforce tail to aid flow control and facilitate easy cleaning
- Tactile feedback upon closing



Nova 1 FoldUp Soft Mouldable Pouch

Box of 10



MIDI	
OPAQUE	STARTER HOLE
881-25	25 – 35 mm (stretch)
881-30	30 – 40 mm (stretch)
881-35	35 – 45 mm (stretch)

Moderma Flex CeraPlus Pouch – with Tape Border

Box of 10

- **Ceramide** Infused Barrier
- Lock 'n Roll Microseal Closure
- AF300 Filter
- Convexity Depth = 4.16 mm



MIDI		
CLEAR	VIEWING OPTION	STARTER HOLE
53700		15 – 38 mm (CTF)
53701		15 – 51 mm (CTF)

MAXI		
CLEAR	VIEWING OPTION	STARTER HOLE
53900	54800	15 – 38 mm (CTF)
53901	54801	15 – 51 mm (CTF)
	54825	25 mm
	54830	30 mm
	54835	35 mm

NovaLife TRE 1 Open Convex Pouch

Box of 10



NEW

MIDI		
CLEAR	EASIVIEW	STARTER HOLE
3042-37	3041-37	15 – 37 mm (CTF)
	3041-46	15 – 46 mm (CTF)
	3045-25	25 mm (PC)
	3041-30	30 mm (PC)
	3041-35	35 mm (PC)
	3041-40	40 mm (PC)
MAXI		
CLEAR	EASIVIEW	STARTER HOLE
3042-59	3041-59	15 – 59 mm (CTF)

NovaLife GX 1 Open Convex Pouch

Box of 10



MINI		
CLEAR	OPAQUE	STARTER HOLE
	941-25	25 mm (PC)
MIDI		
CLEAR	OPAQUE	STARTER HOLE
942-37	941-37	10 – 37 mm (CTF)
942-46	941-46	10 – 46 mm (CTF)
942-30	941-30	30 mm (PC)
942-35	941-35	35 mm (PC)
MAXI		
CLEAR	OPAQUE	STARTER HOLE
942-59	941-59	10 – 59 mm (CTF)
	941-50	50 mm (PC)

Moderma Flex FlexWear Pouch – with Tape Border

Box of 10

- **FlexWear** Durable Skin Barrier
- Lock 'n Roll Microseal Closure
- Belt Tabs
- AF300 Filter
- Convexity Depth = 4.16 mm



MAXI		
CLEAR	OPAQUE	STARTER HOLE
56900	56800	15 – 38 mm (CTF)
56901	56801	15 – 51 mm (CTF)
56920		20 mm
	56825	25 mm
56930	56830	30 mm
56935	56835	35 mm
	56840	40 mm

Moderma Flex SoftFlex Pouch

Box of 30

- **SoftFlex** Skin Barrier
- Lock 'n Roll Microseal Closure
- Belt Tabs
- AF300 Filter
- Convexity Depth = 4.16 mm



MAXI		
CLEAR	OPAQUE	STARTER HOLE
26900	26800	15 – 38 mm (CTF)
	26801	15 – 51 mm (CTF)
26925	26825	25 mm
26930	26830	30 mm
26935	26835	35 mm

NovaLife GX+ 1 Open Convex Pouch

Box of 10



MIDI		
CLEAR	EASIVIEW	STARTER HOLE
	9510-37	10 – 37 mm (CTF)
	9510-25	25 mm (PC)
	9510-30	30 mm (PC)

MAXI		
CLEAR	EASIVIEW	STARTER HOLE
	9510-59	10 – 59 mm (CTF)

Nova 1 FoldUp Convex Pouch

Box of 10



MIDI		
CLEAR	OPAQUE	STARTER HOLE
842-37	841-37	15 – 37 mm (CTF)
842-46	841-46	15 – 46 mm (CTF)
	841-20	20 mm (PC)
842-25	841-25	25 mm (PC)
842-30	841-30	30 mm (PC)
842-35	841-35	35 mm (PC)

Moderma Flex CeraPlus Multi-Chamber Pouch – with Tape Border

Box of 10

- **Ceramide** Infused Barrier
- Enhanced Flow Control Tap
- Multi Chamber Design



NEW
Ceramide
Infused

MAXI	
CLEAR	STARTER HOLE
59400	15 – 64 mm (CTF)

Moderma Flex Multi-Chamber Pouch – with Tape Border

Box of 10

- **Flextend** Skin Barrier
- Enhanced Control Flow Tap
- Multi Chamber Design



MAXI		
TRANSPARENT	OPAQUE	STARTER HOLE
2949064	2909064	13 – 64 mm (CTF)
	2909125	25 mm
	2909130	30 mm
	2909135	35 mm

Moderma Flex Pouch – with Tape Border

Box of 10

- **Flextend** Skin Barrier
- Control Flow Tap



MAXI		
CLEAR	OPAQUE	STARTER HOLE
29400		13 – 64 mm (CTF)
29425		25 mm
29430		30 mm

The novaLife TRE barrier is designed to provide a secure, flexible seal to help protect the skin from stoma output while being easy to remove and gentle to the skin.

- Tough on stoma output, gentle on skin
- Designed for skin protection
- Can help your patients feel comfortable in their own skin

pH balance 

Absorption 

Adhesion 

Designed to help maintain the pH balance of naturally healthy skin

Easy to remove

- TRE barrier is designed to be easy to remove and gentle on the skin

36% Thinner design

- TRE barrier is only 0.7 mm thick compared to our traditional 1.1 mm hydrocolloid

Acid mantle protection

- Maintaining the acid mantle of the skin is essential to help protect its antimicrobial properties

Centred starter hole

- Now with centred starter hole

Flexibility

- Designed to be flexible to allow for natural body movements

Absorption

- Designed to help absorb excess moisture while still maintaining physical integrity

Adhesion

- Designed to have a high tack surface to provide secure adhesion

Designed to absorb stoma output and perspiration

Tapering

- The barrier tapers to 0.3 mm

Friendly to skin

- Unfriendly to digestive enzymes

Two layer barrier designed to prevent curling

Primary layer

- High absorption close to the stoma where protection is needed the most

Secondary layer

- Utilises standard hydrocolloid technology and acts as moisture protection at the edge of the primary layer

[Contents](#)

[Index](#)

Moderna Flex CeraPlus Multi-Chamber Pouch – with Tape Border

Box of 10

- **Ceramide** Infused Barrier
- Enhanced Flow Control Tap
- Multi-Chamber Construction
- With Belt Tabs
- Convexity Depth = 6 mm



NEW
Ceramide
Infused

MAXI		
CLEAR	VIEWING OPTION	STARTER HOLE
59900		15 – 25 mm (CTF)
59901		15 – 38 mm (CTF)
59902		15 – 55 mm (CTF)
	59820	20 mm
	59825	25 mm
	59830	30 mm

Nova 1 Urostomy X3 Pouch

Box of 10



MIDI		
CLEAR	OPAQUE	STARTER HOLE
	891-12	12 – 46 mm (CTF)
	891-20	20 mm (PC)
	891-25	25 mm (PC)
	891-30	30 mm (PC)

Moderma Flex CeraPlus Multi-Chamber Pouch – with Tape Border

Box of 10

- **Ceramide** Infused Barrier
- Enhanced Flow Control Tap
- With Belt Tabs
- Convexity Depth = 4.16 mm



NEW
Ceramide
Infused

MAXI		
CLEAR	VIEWING OPTION	STARTER HOLE
59700		15 – 25 mm (CTF)
59701		15 – 38 mm (CTF)
59702		15 – 51 mm (CTF)
	59620	20 mm
	59625	25 mm
	59630	30 mm

Moderma Flex Multi-Chamber Pouch – with Tape Border

Box of 10

- **Flextend** Skin Barrier
- Enhanced Flow Control Tap
- Convexity Depth = 4.16 mm



IMPROVED
DESIGN

MAXI		
TRANSPARENT	OPAQUE	STARTER HOLE
2979025	2969025	13 – 25 mm (CTF)
2979038	2969038	13 – 38 mm (CTF)
2979120	2969120	20 mm
2979125	2969125	25 mm
2979130	2969130	30 mm
2979135	2969135	35 mm
2979140	2969140	40 mm

Moderma Flex Pouch – with Tape Border

Box of 10

- **Flextend** Skin Barrier
- Belt Tabs
- Control Flow Tap
- Convexity Depth = 4.16 mm



MAXI		
CLEAR	OPAQUE	STARTER HOLE
29700		13 – 25 mm (CTF)
29701		13 – 38 mm (CTF)
29720		20 mm
29725		25 mm
29730		30 mm
29735		35 mm

Nova 1 Urostomy Convex Pouch

Box of 10



MIDI	
CLEAR	STARTER HOLE
894-37	15 – 37 mm (CTF)
894-46	15 – 46 mm (CTF)
894-25	25 mm (PC)
894-30	30 mm (PC)

Non-Sterile Urostomy Night Drainage Bag

Box of 10

- T-Tap drainage outlet
- Connects directly to Hollister and Dansac Urostomy pouches, no need for connectors
- Punch holes for a night stand attachment
- 1200 mm wide bore inlet tubing

CODE	CAPACITY
5550	2000 ml





TWO-PIECE BARRIERS



The CeraPlus Skin Barrier *with Remois Technology**

The Ceramide Advantage

Ceramide is a natural component of human skin, consisting of waxy lipid molecules. Ceramide links cells of the outer skin (epidermis) together to form a waterproof, protective barrier. Ceramides help prevent water loss that can lead to skin damage and dryness.

The CeraPlus Skin Barrier is Designed to Support Healthy Skin From the Start

CeraPlus is designed to support healthy skin from the start. It features a formulation infused with ceramide to help:

- Protect the skin's natural moisture barrier
- Maintain good peristomal skin health
- Minimise moisture loss from the skin (transepidermal water loss – TEWL).

Introducing the NEW Two-Piece CeraPlus Soft Convex Skin Barrier



* Remois is a technology of Alcare Co., Ltd.



Dansac NovaLife TRE

The best skin is healthy skin

NovaLife TRE ostomy skin barriers





Takes skin protection to another level. Three to be exact.

- Available in flat, soft and firm convex barriers
- Stays in place, yet is easy to remove
- Designed to absorb stoma output and perspiration
- Maintaining normal skin pH helps protect it from damaging stoma output
- Two-piece wafers available late 2020



Hollister Colour Matching System

To help you select compatible components, every box label is colour-coded. First, select the proper skin barrier; then select a matching colour-coded pouch.




 44 mm	 70 mm
 57 mm	 102 mm

Two-Piece Barriers – Flat

New Image CeraPlus Barrier

Box of 5

- **Ceramide** Infused Skin Barrier – to help maintain healthy peristomal skin.
- Floating Flange
- CenterPointLock System

TAPE	RING SIZE	STARTER HOLE
11202	 44 mm	12 – 32 mm (CTF)
11203	 57 mm	12 – 44 mm (CTF)
11204	 70 mm	12 – 57 mm (CTF)














**Ceramide
Infused**

New Image Barrier

Box of 5

- **FlexWear** Skin Barrier – Flat
- Floating Flange
- CenterPointLock System

TAPE	NO TAPE	RING SIZE	STARTER HOLE
14202	15202	 44 mm	12 – 32 mm (CTF)
14203	15203	 57 mm	12 – 44 mm (CTF)
14204	15204	 70 mm	12 – 57 mm (CTF)
14206		 102 mm	0 – 89 mm (CTF)
14302		 44 mm	19 mm
14303		 44 mm	22 mm
14304	16404	 44 mm	25 mm
14305		 44 mm	29 mm
14306	16406	 44 mm	32 mm
14307	16407	 57 mm	35 mm
14308		 57 mm	38 mm







Tape Border



No Tape

Dansac Colour Matching System

To help you select compatible components, every box label is colour-coded. First, select the proper skin barrier; then select a matching colour-coded pouch. All Nova and NovaLife 2-Piece products can be used interchangeably.

	36 mm		55 mm
	43 mm		70 mm





Two-Piece Barriers – Flat

NovaLife TRE 2 Flat Wafer

Box of 5







Available late 2020

CODE	RING SIZE	STARTER HOLE
3143-15	 43 mm	15 – 35 mm (CTF)
3155-15	 55 mm	15 – 47 mm (CTF)
3143-25	 43 mm	25 mm (PC)
3143-30	 43 mm	30 mm (PC)

NovaLife GX 2 Flat Wafer

Box of 5













CODE	RING SIZE	STARTER HOLE
1770-10	 70 mm	15 – 62 mm (CTF)
1743-25	 43 mm	25 mm (PC)
1755-35	 55 mm	35 mm (PC)
1770-50	 70 mm	50 mm (PC)

Nova 2 Flat Wafer

Box of 5



CODE	RING SIZE	STARTER HOLE
1143-15	 43 mm	15 – 35 mm (CTF)
1155-15	 55 mm	15 – 47 mm (CTF)
1170-15	 70 mm	15 – 62 mm (CTF)
1136-20	 36 mm	20 mm (PC)
1143-25	 43 mm	25 mm (PC)
1143-30	 43 mm	30 mm (PC)
1155-35	 55 mm	35 mm (PC)
1155-40	 55 mm	40 mm (PC)
1170-45	 70 mm	45 mm (PC)
1170-50	 70 mm	50 mm (PC)

New Image FormaFlex Skin Barrier & Shape-to-Fit – with Tape Border

Box of 5

- Extended Wear Skin Barrier
- No measuring, no tracing and no cutting – stretch to fit
- Floating Flange

CODE	RING SIZE	STARTER HOLE
14102	44 mm	Up to 32 mm (stretch)
14103	57 mm	Up to 43 mm (stretch)
14104	70 mm	Up to 57 mm (stretch)



Nova 2 X3 Wafer

Box of 5



CODE	RING SIZE	STARTER HOLE
1343-15	■ 43 mm	15 – 30 mm (CTF)
1355-15	■ 55 mm	15 – 42 mm (CTF)
1343-25	■ 43 mm	25 mm (PC)
1355-35	■ 55 mm	35 mm (PC)

Nova 2 Soft Mouldable Wafer

Box of 5



CODE	RING SIZE	STARTER HOLE
1443-21	■ 43 mm	21 – 31 mm (stretch)
1443-25	■ 43 mm	25 – 35 mm (stretch)
1443-28	■ 43 mm	28 – 38 mm (stretch)
1455-32	■ 55 mm	32 – 42 mm (stretch)
1455-38	■ 55 mm	38 – 48 mm (stretch)

New Image CeraPlus Barrier – with Tape Border

Box of 5

- **Ceramide** Infused Skin Barrier – to help maintain healthy peristomal skin
- Floating Flange
- Convexity Depth = 6 mm



NEW
Ceramide
Infused

CODE	RING SIZE	STARTER HOLE
11702	44 mm	15 – 25 mm (CTF)
11703	57 mm	15 – 38 mm (CTF)
11704	70 mm	15 – 51 mm (CTF)
11903	44 mm	22 mm
11904	44 mm	25 mm
11905	57 mm	29 mm
11906	57 mm	32 mm

Two-Piece Barriers – Firm Convex

New Image CeraPlus Barrier – with Tape Border

Box of 5

- **Ceramide** Infused Barrier
- Floating Flange
- CenterPointLock System
- Convexity Depth = 4.16 mm



Ceramide
Infused

CODE	RING SIZE	STARTER HOLE
11402	44 mm	15 – 25 mm (CTF)
11403	57 mm	15 – 38 mm (CTF)
11404	70 mm	15 – 51 mm (CTF)
11503	44 mm	22 mm
11504	44 mm	25 mm
11505	57 mm	29 mm
11506	57 mm	32 mm

New Image Barrier

Box of 5

- **FlexWear** Skin Barrier
- Floating Flange
- CenterPointLock System
- Convexity Depth = 4.16 mm



No Tape

TAPE	NO TAPE	RING SIZE	STARTER HOLE
	15403	57 mm	12 – 38 mm (CTF)
	15404	70 mm	15 – 51 mm (CTF)
	15503	44 mm	22 mm
	15504	44 mm	25 mm
	15505	57 mm	29 mm
14907	15506	57 mm	35 mm
14908		57 mm	38 mm

NovaLife TRE 2 Convex Wafer

Box of 5



NEW

Available late 2020

CODE	RING SIZE	STARTER HOLE
3243-15	43 mm	15 – 30 mm (CTF)
3255-15	55 mm	15 – 42 mm (CTF)
3243-25	43 mm	25 mm (PC)
3243-30	43 mm	30 mm (PC)
3255-35	55 mm	35 mm (PC)
3255-40	55 mm	40 mm (PC)

NovaLife GX 2 Convex Wafer

Box of 5



CODE	RING SIZE	STARTER HOLE
1943-10	43 mm	10 – 30 mm (CTF)
1955-10	55 mm	10 – 42 mm (CTF)
1970-10	70 mm	10 – 59 mm (CTF)
1936-20	36 mm	20 mm (PC)
1936-23	36 mm	23 mm (PC)
1943-25	43 mm	25 mm (PC)
1943-30	43 mm	30 mm (PC)
1955-35	55 mm	35 mm (PC)
1955-40	55 mm	40 mm (PC)
1970-45	70 mm	45 mm (PC)

Nova 2 Convex Wafer

Box of 5



CODE	RING SIZE	STARTER HOLE
1543-15	43 mm	15 – 30 mm (CTF)
1555-15	55 mm	15 – 42 mm (CTF)
1543-25	43 mm	25 mm (PC)
1543-30	43 mm	30 mm (PC)
1555-35	55 mm	35 mm (PC)
1555-40	55 mm	40 mm (PC)

New Image FlexWear Barrier – with Tape Border

Box of 5

- **FlexWear** Skin Barrier
- CenterPointLock System
- Floating Flange
- Convexity Depth = 4.16 mm



CODE	RING SIZE	STARTER HOLE
13402	44 mm	15 – 25 mm (CTF)
13403	57 mm	15 – 38 mm (CTF)
13404	70 mm	15 – 51 mm (CTF)
13503	44 mm	22 mm
13504	44 mm	25 mm
13505	57 mm	29 mm
13506	57 mm	32 mm



TWO-PIECE POUCHES

Liberty
freedom to live

A Great Pouch Becomes Even Better

Regardless of what the day brings, Hollister's urostomy pouch helps maximise discretion and comfort.

The improved New Image urostomy pouch from Hollister comes with a multi-chamber design and improved soft touch control flow tap.

- Odour-barrier pouch film
- Transparent and beige options
- Drain valve teardrop indicator
- Belt tabs
- Non-return valves



 Hollister

NovaLife

For confidence and security

Today, people demand functional, aesthetically designed products that suit their lifestyle. We believe that stoma pouches should be no different. Next to discretion, security has the highest priority for people who wear stoma pouches. With NovaLife pouches, we have pushed the limits of design to make them as discreet and secure as possible.

- Low headspace
- Release tab – for simple and secure removal
- Oval shaped barrier – to suit the natural curves of the body
- Highly effective NovaLife filter
- EasiView™ viewing option – allows easy inspection of the stoma and facilitates positioning of the pouch
- Available in TRE, GX and GX+ barriers



dansac 

Hollister Colour Matching System

To help you select compatible components, every box label is colour-coded. First, select the proper skin barrier; then select a matching colour-coded pouch.

44 mm

70 mm57 mm

Two-Piece Pouches – Closed

New Image Closed Pouch

Box of 30

- Belt Tabs
- AF300 Filter





MIDI	
OPAQUE	RING SIZE
18392/30	44 mm
18393/30	57 mm

MAXI	
OPAQUE	RING SIZE
18372	44 mm
18373	57 mm
18374	70 mm

Dansac Colour Matching System

To help you select compatible components, every box label is colour-coded. First, select the proper skin barrier; then select a matching colour-coded pouch. All Nova and NovaLife 2-Piece products can be used interchangeably.


	36 mm		55 mm
	43 mm		70 mm





Two-Piece Pouches – Closed

NovaLife 2 Closed Pouch

Box of 30




MINI	
OPAQUE	RING SIZE
1303-43	 43 mm





MIDI	
OPAQUE	RING SIZE
1301-36	 36 mm
1301-43	 43 mm
1301-55	 55 mm
1301-70	 70 mm

Nova 2 Closed Pouch

Box of 30



MINI	
OPAQUE	RING SIZE
1203-36	 36 mm

MIDI	
OPAQUE	RING SIZE
1201-36	 36 mm
1201-43	 43 mm
1201-55	 55 mm
1201-70	 70 mm

New Image Drainable Pouches

Box of 20 (MINI)

Box of 10 (MAXI)

- Lock 'n Roll Microseal Closure
- Belt Tabs
- AF300 Filter



MINI		
CLEAR	OPAQUE	RING SIZE
	18282	■ 44 mm
	18283	■ 57 mm

MAXI		
CLEAR	OPAQUE	RING SIZE
18192	18182	■ 44 mm
18193	18183	■ 57 mm
18194	18184	■ 70 mm
18006*		■ 102 mm

* No Filter

New Image Drainable Pouches

Box of 10

- Clamp Closure
- Belt Tabs
- AF300 Filter








MAXI		
CLEAR	OPAQUE	RING SIZE
	18142	■ 44 mm
	18143	■ 57 mm
	18144	■ 70 mm

NovaLife 2 Open Pouch

Box of 30






MIDI		
CLEAR	OPAQUE	RING SIZE
	1315-43	 43 mm
	1315-55	 55 mm
	1315-70	 70 mm

MAXI		
CLEAR	OPAQUE	RING SIZE
	1321-43	 43 mm
	1321-55	 55 mm

Nova 2 FoldUp Pouch

Box of 10



MIDI		
CLEAR	OPAQUE	RING SIZE
	1215-43	 43 mm
	1215-55	 55 mm
	1215-70	 70 mm

Two-Piece Pouches – High Output



New Image High Output Drainable Pouch

Box of 10

- No Filter
- Soft Tap
- Belt Tabs



MAXI		
CLEAR	OPAQUE	RING SIZE
18012		44 mm
18013		57 mm
18014		70 mm

Two-Piece Pouches – Urostomy

New Image Urostomy Pouch

Box of 10

- Multi Chamber
- New Control Flow Tap
- Belt Tabs



MAXI		
TRANSPARENT	OPAQUE	RING SIZE
18902	18912	44 mm
18903	18913	57 mm
18904	18914	70 mm

Urostomy Pouch

Box of 10

- Belt Tabs
- Control Flow Tap



MAXI		
TRANSPARENT	OPAQUE	RING SIZE
18402	18412	44 mm
18403	18413	57 mm
18404		70 mm

Non-Sterile Urostomy Night Drainage Bag

Box of 10





CODE	CAPACITY
5550	2000 ml

Nova 2 Urostomy Pouch

Box of 10



MIDI		
CLEAR	OPAQUE	RING SIZE
1218-43	1217-43	 43 mm
1218-55	1217-55	 55 mm



PAEDIATRIC POUCHES

Liberty
freedom to live

PAEDIATRIC
POUCHES

[Contents](#)

[Index](#)

Newborn Pouch/Fistula Pouch – Flat

Box of 15

- **SoftFlex** Skin Barrier
- Use with or without the drain valve
- Soft wire closure
- No filter



ULTRA-CLEAR	BARRIER SIZE	POUCH LENGTH
3778	35 mm x 22 mm (CTF)*	16.5 cm

One-Piece Urostomy – Flat

Box of 10

- **SoftFlex** Skin Barrier



CLEAR	STARTER HOLE	POUCH LENGTH
3797	0 – 38 mm (CTF)*	22 cm

Paediatric Pouches

dansac 

Nova 1 Infant / Paediatric

Box of 30



CLEAR	STARTER HOLE
818-00	0 – 40 mm (CTF)*
818-10	10 – 40 mm (CTF)

* No starter hole



ACCESSORIES

Liberty
freedom to live

[Contents](#)

[Index](#)

Infused with the Science of Skin Health

Presenting a full range of options to help support the healthy peristomal skin your patients deserve

Adapt CeraRing™ barrier rings – delivering the ultimate combination of fit and formulation.

Hollister Ostomy Care is dedicated to delivering products that provide an optimum fit to prevent leakage, along with the right skin barrier formulation for protection and maintenance of healthy skin. This full line of Adapt CeraRing barrier rings offers the benefits of a ceramide-infused formulation and custom-fit shape options to enhance the performance of both one and two-piece pouching systems.



Sealing in Skin Health

Dansac TRE seal is more than just a seal

– with three levels of protection, the Dansac TRE seal has been designed to help keep healthy skin naturally healthy.



Adhesion

Designed to provide a secure, flexible seal to protect the skin from stoma fluid and to be easy to remove.



Absorption

Helps absorb excess moisture to maintain skin's natural balance without the seal losing internal or external strength.



pH balance

Designed to help manage the skin-damaging effects of digestive enzyme activity.



Adapt CeraRing – Standard and Slim

Box of 10

- Can be stretched without flaking or breaking
- Moulded to fit
- Once it has been shaped it keeps its shape



CODE	OUTER DIAMETER	THICKNESS
8805	48 mm	Standard
8815	48 mm	Slim

Adapt CeraRing – Convex

Box of 10



CODE	INNER DIAMETER
89520	20 mm - can be stretched to 25 mm
89530	30 mm - can be stretched to 35 mm
89540	40 mm - can be stretched to 45 mm

Adapt CeraRing – Convex Oval

Box of 10



CODE	INNER DIAMETER
89601	22 x 38 mm
89602	30 x 48 mm
89603	38 x 56 mm

Adapt Flextend M Barrier Ring

Box of 10

- 98 mm for a wider fit all around
- Flat, circular-shaped hydrocolloid skin barriers with a centre opening. Fill uneven skin surfaces around the stoma. Easily shaped to the dimensions of the stoma.



CODE	OUTER DIAMETER
7806	98 mm

Dansac TRE – Seals

Box of 10 (18 mm)

Box of 30 (30 mm & 40 mm)



NEW

CODE	STARTER HOLE	OUTER DIAMETER
072-48	18 mm	48 mm
071-30	30 mm	57.5 mm
071-40	40 mm	70 mm

Dansac Seals

Box of 30



CODE	SIZE
070-20	20 mm
070-30	30 mm
070-40	40 mm

Dansac GX-tra Seals

Box of 30



CODE	SIZE
725-20/30	20 mm
725-30/30	30 mm
725-40/30	40 mm
725-50/30	50 mm

Dansac Soft Paste

Each



CODE	DESCRIPTION
77550-0	50 g tube

Dansac Skin Lotion

Each



CODE	DESCRIPTION
70000-0000	200 ml bottle

Adapt Skin Barrier Paste

Each



CODE	DESCRIPTION	QUANTITY
79300	57 gm tube	Box of 1

Flextend Skin Barrier

Box of 5



CODE	SIZE
8800	10 cm x 10 cm

Adapt No-Sting Skin Protective Wipes

Box of 50 Wipes

- Alcohol free
- Does not 'pill' when drying



CODE	DESCRIPTION
7917	Box 50 Wipes

Adapt No-Sting Adhesive Remover Wipes

Box of 50 Wipes

- For Adhesives and Barriers
- Alcohol free



CODE	DESCRIPTION
7760	Box 50 Wipes

Adapt Stoma Powder

Each



CODE	DESCRIPTION
7906	28.3 g puff bottle

Adapt Skin Conditioning Creme

Each



CODE	DESCRIPTION
7220	120 ml bottle

Dansac Skin Lotion Wipes

Box of 50



CODE	DESCRIPTION
71000-0000	Individual wipes

Dansac EasiSpray™ Adhesive Remover

Each



CODE	DESCRIPTION
083-01	50 ml bottle

Belts

Each



CODE	SIZE
090-02	67 – 100 cm Medium
090-03	113 – 150 cm Large

Dansac Nodor S

Each



CODE	SIZE
080-00	50 ml

Windless Gas Suppressant

Each



CODE	SIZE
W9759	62 Tablets
W9760	250 ml Liquid

Adapt Medical Adhesive Spray

Each

- Medical Adhesive
- Improves adhesive contact between skin and barrier



CODE	DESCRIPTION
7730	90 g spray can

m9 Odour Eliminator Drops

Each



CODE	DESCRIPTION
7715	30 ml bottle

Odour Eliminators and Lubricants

Each



CODE	DESCRIPTION	QUANTITY
78500	236 ml bottle	Each
78501	8 ml sachet	50 Sachet/Box

Ostomy Belts

Each



CODE	SIZE
7300	58 – 109 cm Medium
7299	86 – 165 cm Large

Urostomy Drain Tube Adapter

Box of 10



CODE
7331

Drainable Pouch Clamp

Box of 20



CODE	COLOUR
8770	Beige



Stoma
Skin
Soul



Feeling comfortable
in your own skin

Nils, Crohn's Rebel

Dansac TRE seal is more than just a seal

With three levels of protection, the Dansac TRE seal has been designed to help keep healthy skin naturally healthy.

- Formulated to provide a secure, flexible seal to protect the skin from stoma fluid, while maintaining the pH balance of naturally healthy skin
- Designed to help make everyday life comfortable and deliver the confidence and reassurance you deserve
- Available in 3 sizes to cover a range of stomas




Really, it's just about your patients feeling comfortable in their own skin!

For more information, please call Customer Service
on 0800 678 669 or visit www.dansac.co.nz

Stoma Cone Irrigation Starter Set

- Stoma irrigation drainable
- 12 Stoma irrigation drainable sleeves
- Stoma lubricant 118 ml bottle
- Medium ostomy belt
- Drainable pouch clamp

CODE		
7719		



Cone Irrigator Kit


- Box of 1
- (1) Irrigator
 - (1) Stoma Cone with Connector
 - Cleaning Brush

CODE		
7718		



Stoma Cone with Connector

Each



CODE		
7723		

Irrigator Sleeve with Belt Tabs

Box of 20



CODE	OPENING SIZE	LENGTH
7724	76 mm	56 cm

Dansac Irrigation Set



CODE	DESCRIPTION
900-35	Complete set Hydrocolloid one-piece system

Dansac Irrigation Set Accessories

CODE	DESCRIPTION
950-20	Sleeve with Silicone Ring two-piece system (ea)
950-35	Sleeve Hydrocolloid one-piece system (ea)
95200-0000	Water Bag
95225-0000	Connecting Tube with Regulation Clamp
95205-0000	Cone with Tubing



INTERMITTENT CATHETERS & DRAINAGE BAGS

Liberty
freedom to live

Non-Lubricated Intermittent Catheter

Box of 30

- Smooth eyelets
- Designed for comfort
- **Not** made with natural rubber latex



Tiemann Tip

CODE	SIZE	DESCRIPTION	LENGTH
10616	6Fr	Firm	40 cm
10816	8Fr	Firm	40 cm
11016	10Fr	Firm	40 cm
11216	12Fr	Firm	40 cm
11416	14Fr	Firm	40 cm
11616	16Fr	Firm	40 cm
11816	18Fr	Firm	40 cm
11006	10Fr	Firm	15 cm
11206	12Fr	Firm	15 cm
11406	14Fr	Firm	15 cm
11226	12Fr	Firm Tiemann Tip	40 cm
11426	14Fr	Firm Tiemann Tip	40 cm

Wound Drainage & Night Drainage Bag

Wound Drainage Collectors (Sterile)

Each



CODE	DESCRIPTION	QUANTITY
9701	For Wounds up to 7.6 cm	Each

Non-Sterile Urostomy Night Drainage Bag

Box of 10




CODE	CAPACITY
5550	2000 ml

TUBE ATTACHMENT DEVICES (TAD)

Liberty
freedom to live

The AnchorFast Oral endotracheal tube fastener provides a convenient means to hold an oral endotracheal (ET) tube securely in place without the use of adhesive tape.



AnchorFast Oral Endotracheal Tube Fastener

Box of 12



CODE	DESCRIPTION
9799	Provides a convenient means to hold an oral endotracheal (ET) tube (5.0 – 10.0 mm) securely in place without the use of adhesive tape and is designed to ease access to the oral cavity to help optimise oral care.

AnchorFast Guard Oral Endotracheal Tube Fastener – with the addition of integrated tube protection.

Box of 12



CODE	DESCRIPTION
9800	Integrated tube protection sleeve helps prevent occlusion and designed to ease access to the oral cavity to help optimise oral care. Unlike separate bite blocks, it's there before it's needed.

AnchorFast SlimFit Oral Endotracheal Tube Fastener

Box of 12



**NEW DESIGN TO
ACCOMMODATE
SMALLER FACES**

CODE	DESCRIPTION
9787	Provides a convenient means to hold an oral endotracheal (ET) tube (5.0 - 8.0 mm) securely in place without the use of adhesive tape and is designed to ease access to the oral cavity to help optimise oral care; with the intent to accommodate smaller faces.

AnchorFast Replacement Neck Band Straps

Box of 4



CODE	DESCRIPTION
9788	Used with AnchorFast and AnchorFast Guard oral endotracheal tube fasteners. Fastens with hook and loop closures.

AnchorFast Oral Endotracheal Tube Fastener

Address intubated patient risks – with confidence.

Every patient dependent on a ventilator for respiration is vulnerable. Ventilator-associated risks can result in increased length of stay and increased hospital costs.¹

<div><div>RISKS</div><div>of intubated patients</div></div> <div>CLINICAL BENEFITS</div> <div>of the AnchorFast oral endotracheal tube fastener</div>

1) Restrepo M, et al. "Economic Burden of Ventilator-Associated Pneumonia Based on Total Resource Utilization." Infection Control and Hospital Epidemiology, 2010; 31(5):509-515.

2) Garcia R, et al. "Reducing ventilator-associated pneumonia through advanced oral-dental care: a 48-month study." American Journal of Critical Care 18.6 (2009): 523-32.

3) Miller, Kenneth, et al. "A Utilization of AnchorFast Oral Endotracheal Tube Fastener to Reduce the Incidence of Lip Ulcers." Presented at The 55th International Respiratory Congress. 5-8 December 2009, San Antonio, Texas USA.

4) Veneman, Wade, et al. "An Alternative Approach to Prevent Lip Pressure Ulcers in the Intubated Patient." Presented at the Symposium on Advanced Wound Care. 17-20 April 2010, Orlando, Florida USA.

Tube attachment devices from Hollister provide an easy-to-use alternative to tape for securing tubes, drains, and catheters. Each device features a skin-friendly FlexWear skin barrier to help protect patient skin and an adjustable clamping mechanism to allow for easy tube positioning and better access for site care.

Vertical Tube Attachment Device (VTAD)

Box of 5

- Secures drains, tubes, or catheters (5-40 Fr) vertical to patient's skin
- Secures tubes at exit site
- Sterile, latex-free



Feeding Tube Attachment Device (FTAD)

Box of 20

- Secures nasogastric/ nasointestinal tubes (5-18 Fr)
- Butterfly-shaped skin barrier conforms to patient's nose
- Unique swivel clamp eliminates pressure on the nares
- Latex-free





RESOURCES



Contents

Identifying PMARSI and PMASD	
Skin Complications	84
Hollister Skin Barrier	
Compositions	86
Dansac TRE pH	
Buffering Technology	87
Convex Barrier Options	88
Using Convexity	90
Where a Soft Convex Skin	
Barrier May Improve Fit	92
Hollister eLearning Centre	94
Dansac eLearning Centre	95
Liberty Medical Academy	96

Identifying PMARSI and PMASD Skin Complications

PMARSI

Peristomal Medical Adhesive Related Skin Injury

- An occurrence in which erythema and/or other manifestation of cutaneous abnormality including but not limited to, Vesicle, bulla, erosion or skin tears persists for 30 minutes or more after removal of the adhesive.¹
- A medical adhesive can be defined as a product used to affix an external component (ie, tape, dressing, catheter, electrode, ostomy pouch, or patch) to the skin.¹

For more information: Refer to the Peristomal Skin Assessment Guide: <http://psag.wocn.org>

Skin Stripping



Adhesive skin injury related to the use and/or removal of the adhesive.

Skin Tear



A skin tear is a traumatic wound caused by mechanical forces, including removal of adhesives. Severity may vary by depth (not extending through the subcutaneous layer).

Tension Skin Blister



Adhesive skin injury related to the use and/or removal of the adhesive.

Ruptured Tension Blisters



Adhesive skin injury related to the use and/or removal of the adhesive.

Folliculitis



Inflammation where the skin becomes irritated, especially in the person who must shave routinely to maintain a seal or with removal of the pouch adhesive causing trauma.² Staph aureus is the most frequent microbial culprit leaving an erythematous and pustular area that resembles candidiasis.²

¹) McNichol L, Lund C, Rosen T, Gray M. Medical Adhesives and Patient Safety: State of the Science. Orthopaedic Nurs. 2013; 32(5): 267-281. ²) Bonnie A, & Beck, DE 2008, *Peristomal dermatology*, Clinical Colon Rectal Surgery, vol.21, no. 1, pp. 41–44. ³) Colwell J, et al . MASD Part 3: Peristomal Moisture–Associated Dermatitis and Periwound Moisture–Associated Dermatitis. J Wound Ostomy Continence Nurs. 2011; 38(5):541-553.

PMASD

Moisture Associated Skin Damage

Moisture Associated Skin Damage (MASD) caused by prolonged exposure to various sources of moisture and their contents including:³

- Urine or Stool
- Mucus
- Perspiration
- Saliva

Peristomal Moisture Associated Skin Damage (PMASD)

- is defined as inflammation and denudation of the skin adjacent to a stoma associated with exposure to effluent such as urine or stool.⁴

For more information: Refer to the Peristomal Skin Assessment Guide: <http://psag.wocn.org>

Irritant Contact Dermatitis



Acute - Areas of skin loss and inflammation due to exposure to stool and/or urine. May be the result of leakage or an incorrect sizing of the ostomy barrier around the stoma.

Irritant Contact Dermatitis



Chronic - Areas of skin overgrowth and areas of skin loss/inflammation due to repeated exposure to stool and/or urine or an incorrect sizing of the ostomy barrier around the stoma.

Maceration



A condition caused by excessive amounts of fluid remaining in contact with the skin or the surface of a wound for extended periods.⁵

Fungal Infection



Fungal infection characterised by a confluent centre and distinct satellite lesions. The environment of the skin under skin barrier may enhance the incidence of peristomal skin infections because a warm, dark, moist environment that may be intermittently soiled by urine or stool, provides an ideal setting for microbial growth.

Allergic Dermatitis



Other - Skin develops an allergic reaction to a foreign substance, the body releases inflammatory chemicals and this results in irritated pruritic skin.

Trauma




Pressure injury may result from forces being exerted on the peristomal skin where opposing forces create tissue damage from higher pressures. May be present as localised injury to the peristomal skin and/or underlying tissues, and are often attributed to belts or incorrect skin barrier fit.

4) Gray M, et al. Peristomal Moisture-Associated Skin Damage in Adults With Fecal Ostomies. J Wound Ostomy Continence Nurs. 2013;00(00):1-11.

5) Cutting, F 2001, *The causes and prevention of maceration of the skin*, Nursing Times, <https://www.nursingtimes.net/clinical-archive/tissue-viability/the-causes-and-prevention-of-maceration-of-the-skin-01-11-2001/w.nursingtimes.net/clinical-archive/tissue-viability/the-causes-and-prevention-of-maceration-of-the-skin-01-11-2001/>, Accessed September 2019.

Hollister Skin Barrier Compositions



Mature cells like bricks
Ceramide functions like mortar

CeraPlus Ceramide Infused Skin Barrier


- Helps skin stay healthy from the start:
 - Retains moisture – decrease in Transepidermal Water Loss (TEWL)
 - Replenishes skin's natural barrier



Cotton fibres
Hydrocolloid
Elastomer

SoftFlex Standard Wear Skin Barrier


- Thin and flexible
- Gentle for frequent removal



Hydrocolloid
Cross-linking
Elastomer

FlexWear Extended Wear Skin Barrier

- Cross linkers offer erosion resistance
- Excellent flexibility
- High initial tack
- Durable yet gentle on the skin




Hydrocolloid
Cross-linking
Elastomer

Flexlend Extended Wear Skin Barrier

- Strong adhesion and flexibility
- Highly resistant to erosion
- Swells to fit around stoma


Dansac TRE pH Buffering Technology

Proprietary formulation is designed to maintain the slightly acidic pH level of healthy skin, supporting an environment that is friendly for the skin – and unfriendly to the digestive enzymes, helping stop them in their tracks.




Skin Anatomy:

- Skin's function is to keep moisture in and toxins out
- The skin's acid mantle has a pH between 4.0-6.0
- Peristomal skin health is maintained when skin remains unbroken and the acid mantle intact



Physiology:

- Digestive enzymes in stoma output can be highly corrosive to skin cells
- Enzymes in their optimal pH environment break down protein-rich skin
- The result is Irritant Contact Dermatitis



Goals:

- To help manage the damaging effects of active digestive enzymes on skin
- Ensure a good seal around the stoma
- Select a seal that may reduce enzyme activity by shielding the acid mantle

Dansac NovaLife TRE ostomy barriers and TRE seal are specifically designed to help maintain the pH balance of naturally healthy skin.

Convex Barrier Options

Convex barrier products have been utilised to help provide a secure and predictable seal. Product selection is made by matching stomal protrusion and peristomal contours to the appropriate barriers. The outward curving of the convex skin barrier is designed to interface with the immediate peristomal skin.

This interface depends not only on the depth, rigidity and gradient of the convex skin barrier, but also on the shape and rigidity of the abdominal wall. For example, a flexible convex barrier can bend and move with the body; a firm convex skin barrier does not bend, rather it remains rigid. Appropriate barrier choice is left to the clinician's assessment of the patient and knowledge of ostomy products.

- Selecting the right convex product may require regular reassessment of products used, condition of the peristomal skin, and the pouching system wear time
- Some ostomy care nurses recommend an ostomy belt to enhance the effect of the convexity.

Barrier Example

Possible Clinical Case Scenarios



Flat Skin Barrier



Flat Skin Barrier

- Flat peristomal plane
- Adequate stoma protrusion
- Centred stoma opening



Flat Skin Barrier with Round or Oval Convex Barrier Ring



Flexible Convexity

- Stoma flush with skin
- Uneven peristomal plane
- Firm abdomen



Integrated Soft Convex Barrier



Soft Convexity

- Gentle pressure required to peristomal plane
- Flush stoma
- Peristomal injury



Integrated Convex Barrier



Firm Convexity

- Creases or folds in peristomal plane
- Off-centred stoma opening
- Stoma retracted below skin surface
- Soft abdomen



Integrated Convex Barrier with Round or Oval Convex Barrier Ring



Firm Convexity with More Depth

- Recessed or flush stoma
- Uneven abdominal contours
- Off-centred stoma opening

1) J.C. Colwell, M.T. Goldberg, & J.E. Carmel (Eds). (2015). WOCN Society Core Curriculum Ostomy Management. Chapter 10. Philadelphia, PA: Wolters Kluwer.

2) Hoefloek, J et al. Use of convexity in ostomy care: A Comprehensive Review. JWOCN, January/February 2017.

Using Convexity

Convexity is the outward curving of the skin barrier and is designed to interface with the immediate peristomal skin. Convex skin barriers promote a good fit between the ostomy pouching system and the peristomal skin where flat barriers would be unsuccessful. This is achieved as the convexity barrier pushes on the surrounding skin, opens or flattens skin folds, fills uneven surfaces or helps the stoma protrude more.



Stoma that is below skin level



Convex skin barrier in place

General Information

- A secure and reliable seal between the peristomal skin and ostomy skin barrier is necessary to protect healthy peristomal skin from stoma discharge
- Skin barriers are either convex or flat. Convex simply describes the shape of the barrier
- Convexity may be recommended to increase the wear time of the ostomy pouching system
- Convexity may be firm or soft

Assessing for Convexity

- A simple patient assessment methodology can assist in determining the need for convexity. Refer to the Hollister **FIT Indicator Tool** (www.hollister.co.nz/education)
- Assessing when convexity is required should be done in consultation with an ostomy care nurse
- Ideally the patient should be assessed in the sitting position
- The pouch should be removed prior to assessment
- A distal stoma assessment is recommended for all loop stomas



Off Centred Stoma Opening



Retracted Stoma



Flush Stoma



Uneven/Irregular



Flaccid Abdomen

Convexity May be Considered When:

- Stoma Assessment
 - Liquid/Loose stoma output
 - The opening of the stoma is at or below the skin surface
 - The opening of the stoma is off-centre
 - The stoma protrudes $<3/4"$ (20 mm) or is retracted below the skin surface
 - The stoma becomes flush/retracted during peristalsis (telescopes)
- Peristomal Assessment
 - The peristomal topography is uneven or irregular
 - The abdominal tone is flaccid



Soft integrated convex barrier



Soft convexity created by adding a convex barrier ring to flat skin barrier

Convex Products

- Convexity can either be built into the pouching system (integrated) or added with an accessory to the skin barrier
- Convexity is available in one-piece and two-piece options
- Convexity is available in pre-sized or cut-to-fit options
 - Note: When using cut-to-fit options for oval or irregular shape stomas, choose a barrier that is closest to the largest measurement so the convex curvature of the barrier is as close to the stoma as possible
- Soft convexity can be integrated into the pouching system or created by adding a convex barrier ring to a flat one-piece pouching system
- Firm convexity can be integrated into the pouching system or created by adding a convex barrier ring to a convex barrier to deepen the curvature when necessary



Firm integrated convex barrier



Firm integrated convex barrier with convex barrier ring to deepen the curvature

Additional Considerations

- The following peristomal conditions do not preclude the use of convex products, however, caution and frequent reassessment is recommended in patients with, Pyoderma Gangrenosum, Caput medusa (peristomal varices), Pressure injury, Crohn's ulcers, Parastomal hernia
- Patients should be encouraged to contact an ostomy care nurse if there are stomal/peristomal changes, pouching system leakage or weight gain/loss

References: • "Patient Assessment Guidelines for Convexity" – Hollister Incorporated © 2013. • "Developing Patient Assessment Guidelines for Convexity" – Drolshagen et al. 2014 Poster Presentation, WOCN Conference, Nashville TN USA. • "Convex Pouching Systems: Best Practice Document for Patients - WOCN 2012. • Use of Convexity in Pouching – A Comprehensive Review – Hoefloek J, Kittscha J, Purnell P. JWOCN 2012:40 (5):506-512. • "Adapt Convex Barrier Rings" – Ostomy Care Tip – Hollister Incorporated © 2012. • Lyon, C. & Smith, A. 2010 – Abdominal Stomas and Their Skin Disorders.

Where a Soft Convex Skin Barrier May Improve Fit



Case 1: Parastomal hernia with flush stoma¹

A parastomal hernia may create a bulge on the abdomen. When convexity is clinically indicated, a firm rigid convex barrier may pose the potential for pressure related skin damage.

A flexible soft convex barrier may be an appropriate solution.



Case 2: Firm abdomen with pressure injury from use of firm rigid convexity^{2,4}

Pressure ulcers may be more likely when a firm convex barrier is pressing against the skin for a prolonged period of time, especially with the addition of an ostomy belt.

A flexible soft convex barrier may provide the correct fit while removing the cause of pressure.



Case 3: Stoma located in a crease³

A firm rigid convex barrier may not conform to the abdominal contours and “pop off” when used in a creased area.

A flexible soft convex barrier may be considered a more appropriate fit.



Case 4: Stoma located in abdominal folds^{2,3}

Abdominal folds can compromise the seal of the barrier. A convex shape can enhance the barrier fit.

A flexible soft convex barrier may match to the correct depth of the folds, conform to the abdominal contours and provide less peristomal pressure.

1) Turnbull, G. The Ostomy Files: Parastomal Hernia. *Ostomy Wound Management* Volume 49- Issue 11 – November 2003

2) J.C. Colwell, M.T. Goldberg, & J.E. Carmel (Eds.). (2004). *Fecal & Urinary Diversions: Management Principles*. In Chapter 12 and Chapter 14. St. Louis, MO: Elsevier Mosby

3) J.C. Colwell, M.T. Goldberg, & J.E. Carvel (Eds.). (2015). *WOCN® Society Core Curriculum Ostomy Management*. In Chapter 10. Philadelphia, PA: Wolters Kluwer



Case 5: Stoma height less than 20 mm (2 cm)³

A stoma that does not protrude above the skin may cause leakage problems under the skin barrier.

A flexible soft convex barrier may provide the right depth of convexity to help with stoma protrusion.



Case 6: Immediate post-op stoma with firm distended abdomen and off-centred lumen at risk for mucocutaneous separation^{2,5}

A mucocutaneous separation may occur as a result of poor healing, infection, or excessive tension at the mucocutaneous junction.

A flexible soft convex barrier may help achieve a correct fit with less pressure at the base of the stoma.



Case 7: Loop stoma³

The distal limb (arrow) of a loop stoma may discharge mucous which can undermine the barrier seal.

A flexible soft convex barrier may help provide the right fit with less pressure around the stoma to minimise undermining.



Case 8: Pyoderma gangrenosum^{2,4}

Trauma to the peristomal skin may initiate and aggravate a pyoderma gangrenosum ulcer. Efforts should be made to alleviate pressure and friction.

A flexible soft convex barrier may provide less pressure than firm convexity.

4) J.E. Carmel, J.C. Colwell, M.T. Goldberg (Eds.), *WOCN Society Core Curriculum Ostomy Management* p 181, 186. Philadelphia: Wolters Kluwer. 2016
 5) External Stoma and Peristomal Complications following Radical Cystectomy and Ileal Conduit Diversion: A Systematic Review - Szymanski, K.M., St-Cyr, D., Alam, T., Kassouf, W. *WOCN Society Core Curriculum Ostomy Management*. p. 192. Philadelphia, PA: Wolters Kluwer. 2010



Looking for
more information
around living with
an ostomy?

Feel confident about equipping your patients with pre and post surgery information to help them live their best life.

Visit the online learning centre – a wealth of ostomy related information and hints and tips on living your best life.

www.hollister.co.nz/ostomylearningcentre



Understanding
an Ostomy



Maintaining
Skin Health



Living With
an Ostomy



Using Ostomy
Products



You can view/download a video on the following topics:

Managing Your Colostomy
Managing Your Ileostomy
Managing Your Urostomy

www.vimeo.com/hollisteranz/colostomy
www.vimeo.com/hollisteranz/ileostomy
www.vimeo.com/hollisteranz/Urostomy



Supporting You and Your Patients with Dansac's 'Living with a Stoma' Online Learning Centre

Having a stoma is a big change in a person's life. How you deal with it depends on your situation, and the information and support you need may differ from what others need.

Visit the online learning centre for educational resources that can help make life with a stoma better – before and after your surgery, and as you get back into your regular routine.

www.dansac.co.nz/livingwithastoma



Facts About Your Stoma

Get information on the different kinds of stomas, and on why someone may need stoma surgery.



Recovery After Stoma Surgery

Learn what to expect after stoma surgery, and how to ease back into your day-to-day life.



Life With Your Stoma

Get tips that can help you take control of your new life and live well with a stoma.



Healthy Skin & Your Stoma

Learn how keeping the skin around your stoma healthy can help enhance the quality of your life.

Visit nurse.libertymed.co.nz to access:



Liberty Medical Academy
Nurse webcasts & online learning



Patient & Nurse Support
Support materials, downloadable literature & customer care



Nurse Portal
News & events for your speciality, product updates & clinical abstracts



Peer2Peer
Interactive peer education programs



Simple Ordering
Custom samples, easy access 24/7



Innovations
Extensive, advanced product portfolios

Better resources, better outcomes.
Supporting your trusted Ostomy Brands.



dansac 



 **Hollister.**



At Liberty Medical, we are committed to fulfilling the many and varied individual needs of people with stomas through an extensive products and supporting materials.

Our patient website contains a range of useful resources for ostomates including the basics, health and wellbeing, hints and tips, personal stories and general information.



www.libertymed.co.nz



NOTES

NOTES

NOTES



INDEX

Liberty
freedom to live

INDEX

[Contents](#)

[Index](#)

Hollister Code Index

Code	Name	Page No.
3184	One-Piece Closed – Flat – Stoma Cap	10
3192	One-Piece Closed – Flat – Compact Stoma Cap	10
3193	One-Piece Closed – Flat – Compact Stoma Cap	10
3778	Pouchkins Paediatric Pouches – Flat SoftFlex Skin Barrier Newborn Pouch/Fistula Pouch	57
3797	Pouchkins Paediatric Pouches– Flat SoftFlex Skin Barrier One-Piece Urostomy	57
5550	Drainage Bags Non-Sterile Urostomy Night Drainage Bag	34
5550	Two-Piece Pouches – Urostomy Non-Sterile Urostomy Night Drainage Bag	52
5550	Wound Drainage & Night Drainage Bag Non-Sterile Urostomy Night Drainage Bag	73
7220	Accessories Adapt Skin Conditioning Crème	64
7299	Accessories Ostomy Belts	66
7300	Accessories Ostomy Belts	66
7331	Accessories Urostomy Drain Tube Adaptor	66
7715	Accessories m9 Odour Eliminator Drops	66
7718	Irrigation Products Cone Irrigator Kit	68
7719	Irrigation Products Stoma Care Irrigation Starter Set	68
7723	Irrigation Products Stoma Cone with Connector	68
7724	Irrigation Products Irrigator Sleeve with Belt Tabs	68
7730	Accessories Adapt Medical Adhesive Spray	66
7760	Accessories Adapt No-Sting Adhesive Remover Wipes	64
7806	Accessories Adapt Flextend M Barrier Ring	62
7906	Accessories Adapt Stoma Powder	64
7917	Accessories Adapt No-Sting Adhesive Protective Wipes	64
8770	Accessories Drainable Pouch Clamp	66
8800	Accessories Flextend Skin Barrier	64
8805	Accessories Adapt CeraRing - Standard and Slim	62
8815	Accessories Adapt CeraRing - Standard and Slim	62
9701	Wound Drainage & Night Drainage Bag Wound Drainage Collectors (Sterile)	73
9782	Tube Attachment Devices (TAD) Vertical Tube Attachment Device (VTAD)	80
9786	Tube Attachment Devices (TAD) Feeding Tube Attachment Device (FTAD)	80
9787	Tube Attachment Devices (TAD) AnchorFast SlimFit Oral Endotracheal Tube Fastener	78
9788	Tube Attachment Devices (TAD) AnchorFast Replacement Neck Band Straps	78
9799	Tube Attachment Devices (TAD) AnchorFast Oral Endotracheal Tube Fastener	78
9800	Tube Attachment Devices (TAD) AnchorFast Guard Oral Endotracheal Tube Fastener – with the addition of integrated tube protection	78
10616	Apogee™ Intermittent Catheter Non-Lubricated Intermittent Catheter	73
10816	Apogee™ Intermittent Catheter Non-Lubricated Intermittent Catheter	73
11006	Apogee™ Intermittent Catheter Non-Lubricated Intermittent Catheter	73
11016	Apogee™ Intermittent Catheter Non-Lubricated Intermittent Catheter	73

Code	Name	Page No.
11202	Two-Piece Barriers – Flat New Image CeraPlus Barrier	38
11203	Two-Piece Barriers – Flat New Image CeraPlus Barrier	38
11204	Two-Piece Barriers – Flat New Image CeraPlus Barrier	38
11206	Apogee™ Intermittent Catheter Non-Lubricated Intermittent Catheter	73
11216	Apogee™ Intermittent Catheter Non-Lubricated Intermittent Catheter	73
11226	Apogee™ Intermittent Catheter Non-Lubricated Intermittent Catheter	73
11402	Two-Piece Barriers – Firm Convex New Image CeraPlus Barrier – with Tape Border	42
11403	Two-Piece Barriers – Firm Convex New Image CeraPlus Barrier – with Tape Border	42
11404	Two-Piece Barriers – Firm Convex New Image CeraPlus Barrier – with Tape Border	42
11406	Apogee™ Intermittent Catheter Non-Lubricated Intermittent Catheter	73
11416	Apogee™ Intermittent Catheter Non-Lubricated Intermittent Catheter	73
11426	Apogee™ Intermittent Catheter Non-Lubricated Intermittent Catheter	73
11503	Two-Piece Barriers – Firm Convex New Image CeraPlus Barrier – with Tape Border	42
11504	Two-Piece Barriers – Firm Convex New Image CeraPlus Barrier – with Tape Border	42
11505	Two-Piece Barriers – Firm Convex New Image CeraPlus Barrier – with Tape Border	42
11506	Two-Piece Barriers – Firm Convex New Image CeraPlus Barrier – with Tape Border	42
11616	Apogee™ Intermittent Catheter Non-Lubricated Intermittent Catheter	73
11702	Two-Piece Barriers – Soft Convex New Image CeraPlus Barrier – with Tape Border	42
11703	Two-Piece Barriers – Soft Convex New Image CeraPlus Barrier – with Tape Border	42
11704	Two-Piece Barriers – Soft Convex New Image CeraPlus Barrier – with Tape Border	42
11816	Apogee™ Intermittent Catheter Non-Lubricated Intermittent Catheter	73
11903	Two-Piece Barriers – Soft Convex New Image CeraPlus Barrier – with Tape Border	42
11904	Two-Piece Barriers – Soft Convex New Image CeraPlus Barrier – with Tape Border	42
11905	Two-Piece Barriers – Soft Convex New Image CeraPlus Barrier – with Tape Border	42
11906	Two-Piece Barriers – Soft Convex New Image CeraPlus Barrier – with Tape Border	42
13402	Two-Piece Barriers – Firm Convex New Image FlexWear Barrier – with Tape Border	44
13403	Two-Piece Barriers – Firm Convex New Image FlexWear Barrier – with Tape Border	44
13404	Two-Piece Barriers – Firm Convex New Image FlexWear Barrier – with Tape Border	44
13503	Two-Piece Barriers – Firm Convex New Image FlexWear Barrier – with Tape Border	44
13504	Two-Piece Barriers – Firm Convex New Image FlexWear Barrier – with Tape Border	44
13505	Two-Piece Barriers – Firm Convex New Image FlexWear Barrier – with Tape Border	44
13506	Two-Piece Barriers – Firm Convex New Image FlexWear Barrier – with Tape Border	44
14102	FormaFlex – Shape to Fit New Image FormaFlex Skin Barrier & Shape-to-Fit – with Tape Border	40
14103	FormaFlex – Shape to Fit New Image FormaFlex Skin Barrier & Shape-to-Fit – with Tape Border	40
14104	FormaFlex – Shape to Fit New Image FormaFlex Skin Barrier & Shape-to-Fit – with Tape Border	40
14202	Two-Piece Barriers – Flat New Image Barrier	38
14203	Two-Piece Barriers – Flat New Image Barrier	38
14204	Two-Piece Barriers – Flat New Image Barrier	38
14206	Two-Piece Barriers – Flat New Image Barrier	38
14302	Two-Piece Barriers – Flat New Image Barrier	38

Code	Name	Page No.
14303	Two-Piece Barriers – Flat New Image Barrier	38
14304	Two-Piece Barriers – Flat New Image Barrier	38
14305	Two-Piece Barriers – Flat New Image Barrier	38
14306	Two-Piece Barriers – Flat New Image Barrier	38
14307	Two-Piece Barriers – Flat New Image Barrier	38
14308	Two-Piece Barriers – Flat New Image Barrier	38
14907	Two-Piece Barriers – Firm Convex New Image Barrier	42
14908	Two-Piece Barriers – Firm Convex New Image Barrier	42
15202	Two-Piece Barriers – Flat New Image Barrier	38
15203	Two-Piece Barriers – Flat New Image Barrier	38
15204	Two-Piece Barriers – Flat New Image Barrier	38
15403	Two-Piece Barriers – Firm Convex New Image Barrier	42
15404	Two-Piece Barriers – Firm Convex New Image Barrier	42
15503	Two-Piece Barriers – Firm Convex New Image Barrier	42
15504	Two-Piece Barriers – Firm Convex New Image Barrier	42
15505	Two-Piece Barriers – Firm Convex New Image Barrier	42
15506	Two-Piece Barriers – Firm Convex New Image Barrier	42
16404	Two-Piece Barriers – Flat New Image Barrier	38
16406	Two-Piece Barriers – Flat New Image Barrier	38
16407	Two-Piece Barriers – Flat New Image Barrier	38
18006	Two-Piece Pouches – Drainable New Image Drainable Pouches	50
18012	Two-Piece Pouches – High Output New Image High Output Drainable Pouch	52
18013	Two-Piece Pouches – High Output New Image High Output Drainable Pouch	52
18014	Two-Piece Pouches – High Output New Image High Output Drainable Pouch	52
18142	Two-Piece Pouches – Drainable New Image Drainable Pouches	50
18143	Two-Piece Pouches – Drainable New Image Drainable Pouches	50
18144	Two-Piece Pouches – Drainable New Image Drainable Pouches	50
18182	Two-Piece Pouches – Drainable New Image Drainable Pouches	50
18183	Two-Piece Pouches – Drainable New Image Drainable Pouches	50
18184	Two-Piece Pouches – Drainable New Image Drainable Pouches	50
18192	Two-Piece Pouches – Drainable New Image Drainable Pouches	50
18193	Two-Piece Pouches – Drainable New Image Drainable Pouches	50
18194	Two-Piece Pouches – Drainable New Image Drainable Pouches	50
18282	Two-Piece Pouches – Drainable New Image Drainable Pouches	50
18283	Two-Piece Pouches – Drainable New Image Drainable Pouches	50
18372	Two-Piece Pouches – Closed New Image Closed Pouch	48
18373	Two-Piece Pouches – Closed New Image Closed Pouch	48
18374	Two-Piece Pouches – Closed New Image Closed Pouch	48
18392/30	Two-Piece Pouches – Closed New Image Closed Pouch	48
18393/30	Two-Piece Pouches – Closed New Image Closed Pouch	48
18402	Two-Piece Pouches – Urostomy Urostomy Pouch	52
18403	Two-Piece Pouches – Urostomy Urostomy Pouch	52

Code	Name	Page No.
18404	Two-Piece Pouches – Urostomy Urostomy Pouch	52
18412	Two-Piece Pouches – Urostomy Urostomy Pouch	52
18413	Two-Piece Pouches – Urostomy Urostomy Pouch	52
18902	Two-Piece Pouches – Urostomy New Image Urostomy Pouch	52
18903	Two-Piece Pouches – Urostomy New Image Urostomy Pouch	52
18904	Two-Piece Pouches – Urostomy New Image Urostomy Pouch	52
18912	Two-Piece Pouches – Urostomy New Image Urostomy Pouch	52
18913	Two-Piece Pouches – Urostomy New Image Urostomy Pouch	52
18914	Two-Piece Pouches – Urostomy New Image Urostomy Pouch	52
22100	One-Piece Closed – Flat Moderma Flex Pouch	8
22120	One-Piece Closed – Flat Moderma Flex Pouch	8
22125	One-Piece Closed – Flat Moderma Flex Pouch	8
22130	One-Piece Closed – Flat Moderma Flex Pouch	8
22135	One-Piece Closed – Flat Moderma Flex Pouch	8
22140	One-Piece Closed – Flat Moderma Flex Pouch	8
22300	One-Piece Closed – Flat Moderma Flex Pouch	8
22301	One-Piece Closed – Flat Moderma Flex Oval Pouch	10
22325	One-Piece Closed – Flat Moderma Flex Pouch	8
22330	One-Piece Closed – Flat Moderma Flex Pouch	8
22335	One-Piece Closed – Flat Moderma Flex Pouch	8
22400	One-Piece Closed – Flat Moderma Flex Pouch	8
22401	One-Piece Closed – Flat Moderma Flex Oval Pouch	10
22500	One-Piece Closed – Flat Moderma Flex Pouch	8
22800	One-Piece Closed – Firm Convex Moderma Flex Pouch	14
22801	One-Piece Closed – Firm Convex Moderma Flex Pouch	14
22820	One-Piece Closed – Firm Convex Moderma Flex Pouch	14
22825	One-Piece Closed – Firm Convex Moderma Flex Pouch	14
22830	One-Piece Closed – Firm Convex Moderma Flex Pouch	14
22835	One-Piece Closed – Firm Convex Moderma Flex Pouch	14
22900	One-Piece Closed – Firm Convex Moderma Flex Pouch	14
22901	One-Piece Closed – Firm Convex Moderma Flex Pouch	14
22925	One-Piece Closed – Firm Convex Moderma Flex Pouch	14
22930	One-Piece Closed – Firm Convex Moderma Flex Pouch	14
22935	One-Piece Closed – Firm Convex Moderma Flex Pouch	14
26064	One-Piece Drainable – Flat Moderma Flex High Output Pouch	18
28300	One-Piece Drainable – Flat Moderma Flex Pouch	16
28301	One-Piece Drainable – Flat Moderma Flex Oval Pouch	18
28325	One-Piece Drainable – Flat Moderma Flex Pouch	16
28330	One-Piece Drainable – Flat Moderma Flex Pouch	16
28335	One-Piece Drainable – Flat Moderma Flex Pouch	16
28340	One-Piece Drainable – Flat Moderma Flex Pouch	16
28400	One-Piece Drainable – Flat Moderma Flex Pouch	16

Code	Name	Page No.
28401	One-Piece Drainable – Flat Moderma Flex Oval Pouch	18
29400	One-Piece Urostomy – Flat Moderma Flex Pouch – with Tape Border	28
29425	One-Piece Urostomy – Flat Moderma Flex Pouch – with Tape Border	28
29430	One-Piece Urostomy – Flat Moderma Flex Pouch – with Tape Border	28
29700	One-Piece Urostomy – Firm Convex Moderma Flex Pouch – with Tape Border	32
29701	One-Piece Urostomy – Firm Convex Moderma Flex Pouch – with Tape Border	32
29720	One-Piece Urostomy – Firm Convex Moderma Flex Pouch – with Tape Border	32
29725	One-Piece Urostomy – Firm Convex Moderma Flex Pouch – with Tape Border	32
29730	One-Piece Urostomy – Firm Convex Moderma Flex Pouch – with Tape Border	32
29735	One-Piece Urostomy – Firm Convex Moderma Flex Pouch – with Tape Border	32
42100	One-Piece Closed – Flat Moderma Flex Quietwear Pouch	8
42125	One-Piece Closed – Flat Moderma Flex Quietwear Pouch	8
42130	One-Piece Closed – Flat Moderma Flex Quietwear Pouch	8
42135	One-Piece Closed – Flat Moderma Flex Quietwear Pouch	8
42140	One-Piece Closed – Flat Moderma Flex Quietwear Pouch	8
42145	One-Piece Closed – Flat Moderma Flex Quietwear Pouch	8
52100	One-Piece Closed – Flat Moderma Flex CeraPlus Pouch	8
52125	One-Piece Closed – Flat Moderma Flex CeraPlus Pouch	8
52130	One-Piece Closed – Flat Moderma Flex CeraPlus Pouch	8
52135	One-Piece Closed – Flat Moderma Flex CeraPlus Pouch	8
52600	One-Piece Closed – Firm Convex Moderma Flex CeraPlus Pouch	14
52601	One-Piece Closed – Firm Convex Moderma Flex CeraPlus Pouch	14
52625	One-Piece Closed – Firm Convex Moderma Flex CeraPlus Pouch	14
52630	One-Piece Closed – Firm Convex Moderma Flex CeraPlus Pouch	14
52635	One-Piece Closed – Firm Convex Moderma Flex CeraPlus Pouch	14
53100	One-Piece Drainable – Soft Convex Moderma Flex CeraPlus Pouch – with Tape Border	20
53101	One-Piece Drainable – Soft Convex Moderma Flex CeraPlus Pouch – with Tape Border	20
53400	One-Piece Drainable – Flat Moderma Flex CeraPlus Pouch – with Tape Border	16
53500	One-Piece Drainable – Flat Moderma Flex CeraPlus Pouch – with Tape Border	16
53700	One-Piece Drainable – Firm Convex Moderma Flex CeraPlus Pouch – with Tape Border	24
53701	One-Piece Drainable – Firm Convex Moderma Flex CeraPlus Pouch – with Tape Border	24
53900	One-Piece Drainable – Firm Convex Moderma Flex CeraPlus Pouch – with Tape Border	24
53901	One-Piece Drainable – Firm Convex Moderma Flex CeraPlus Pouch – with Tape Border	24
54100	One-Piece Drainable – Soft Convex Moderma Flex CeraPlus Pouch – with Tape Border	20
54101	One-Piece Drainable – Soft Convex Moderma Flex CeraPlus Pouch – with Tape Border	20
54125	One-Piece Drainable – Soft Convex Moderma Flex CeraPlus Pouch – with Tape Border	20
54130	One-Piece Drainable – Soft Convex Moderma Flex CeraPlus Pouch – with Tape Border	20
54135	One-Piece Drainable – Soft Convex Moderma Flex CeraPlus Pouch – with Tape Border	20
54200	One-Piece Drainable – Flat Moderma Flex CeraPlus Pouch – with Tape Border	16
54300	One-Piece Drainable – Flat Moderma Flex CeraPlus Pouch – with Tape Border	16
54325	One-Piece Drainable – Flat Moderma Flex CeraPlus Pouch – with Tape Border	16
54330	One-Piece Drainable – Flat Moderma Flex CeraPlus Pouch – with Tape Border	16

Code	Name	Page No.
54335	One-Piece Drainable – Flat Moderma Flex CeraPlus Pouch – with Tape Border	16
54800	One-Piece Drainable – Firm Convex Moderma Flex CeraPlus Pouch – with Tape Border	24
54801	One-Piece Drainable – Firm Convex Moderma Flex CeraPlus Pouch – with Tape Border	24
54825	One-Piece Drainable – Firm Convex Moderma Flex CeraPlus Pouch – with Tape Border	24
54830	One-Piece Drainable – Firm Convex Moderma Flex CeraPlus Pouch – with Tape Border	24
54835	One-Piece Drainable – Firm Convex Moderma Flex CeraPlus Pouch – with Tape Border	24
56800	One-Piece Drainable – Firm Convex Moderma Flex FlexWear Pouch – with Tape Border	26
56801	One-Piece Drainable – Firm Convex Moderma Flex FlexWear Pouch – with Tape Border	26
56825	One-Piece Drainable – Firm Convex Moderma Flex FlexWear Pouch – with Tape Border	26
56830	One-Piece Drainable – Firm Convex Moderma Flex FlexWear Pouch – with Tape Border	26
56835	One-Piece Drainable – Firm Convex Moderma Flex FlexWear Pouch – with Tape Border	26
56840	One-Piece Drainable – Firm Convex Moderma Flex FlexWear Pouch – with Tape Border	26
56900	One-Piece Drainable – Firm Convex Moderma Flex FlexWear Pouch – with Tape Border	26
56901	One-Piece Drainable – Firm Convex Moderma Flex FlexWear Pouch – with Tape Border	26
56920	One-Piece Drainable – Firm Convex Moderma Flex FlexWear Pouch – with Tape Border	26
56930	One-Piece Drainable – Firm Convex Moderma Flex FlexWear Pouch – with Tape Border	26
56935	One-Piece Drainable – Firm Convex Moderma Flex FlexWear Pouch – with Tape Border	26
58300	One-Piece Drainable – Flat Moderma Flex FlexWear Pouch – with Tape Border	16
58325	One-Piece Drainable – Flat Moderma Flex FlexWear Pouch – with Tape Border	16
58400	One-Piece Drainable – Flat Moderma Flex FlexWear Pouch – with Tape Border	16
59400	One-Piece Urostomy – Flat Moderma Flex CeraPlus Multi-Chamber Pouch – with Tape Border	28
59620	One-Piece Urostomy – Firm Convex Moderma Flex CeraPlus Multi-Chamber Pouch – with Tape Border	32
59625	One-Piece Urostomy – Firm Convex Moderma Flex CeraPlus Multi-Chamber Pouch – with Tape Border	32
59630	One-Piece Urostomy – Firm Convex Moderma Flex CeraPlus Multi-Chamber Pouch – with Tape Border	32
59700	One-Piece Urostomy – Firm Convex Moderma Flex CeraPlus Multi-Chamber Pouch – with Tape Border	32
59701	One-Piece Urostomy – Firm Convex Moderma Flex CeraPlus Multi-Chamber Pouch – with Tape Border	32
59702	One-Piece Urostomy – Firm Convex Moderma Flex CeraPlus Multi-Chamber Pouch – with Tape Border	32
59820	One-Piece Urostomy – Soft Convex Moderma Flex CeraPlus Multi-Chamber Pouch – with Tape Border	30
59825	One-Piece Urostomy – Soft Convex Moderma Flex CeraPlus Multi-Chamber Pouch – with Tape Border	30
59830	One-Piece Urostomy – Soft Convex Moderma Flex CeraPlus Multi-Chamber Pouch – with Tape Border	30
59900	One-Piece Urostomy – Soft Convex Moderma Flex CeraPlus Multi-Chamber Pouch – with Tape Border	30
59901	One-Piece Urostomy – Soft Convex Moderma Flex CeraPlus Multi-Chamber Pouch – with Tape Border	30
59902	One-Piece Urostomy – Soft Convex Moderma Flex CeraPlus Multi-Chamber Pouch – with Tape Border	30

Code	Name	Page No.
78500	Accessories Odour Eliminators and Lubricants	66
78501	Accessories Odour Eliminators and Lubricants	66
79300	Accessories Adapt Skin Barrier Paste	64
89004	One-Piece Drainable – Flat Post-Operative Pouch	18
89520	Accessories Adapt CeraRing – Convex	62
89530	Accessories Adapt CeraRing – Convex	62
89540	Accessories Adapt CeraRing – Convex	62
89601	Accessories Adapt CeraRing – Convex – Oval	62
89602	Accessories Adapt CeraRing – Convex – Oval	62
89603	Accessories Adapt CeraRing – Convex – Oval	62
2261038	One-Piece Closed – Soft Convex Moderma Flex Pouch	12
2261055	One-Piece Closed – Soft Convex Moderma Flex Pouch	12
2261125	One-Piece Closed – Soft Convex Moderma Flex Pouch	12
2261130	One-Piece Closed – Soft Convex Moderma Flex Pouch	12
2261135	One-Piece Closed – Soft Convex Moderma Flex Pouch	12
2261140	One-Piece Closed – Soft Convex Moderma Flex Pouch	12
2281025	One-Piece Closed – Soft Convex Moderma Flex Pouch	12
2281038	One-Piece Closed – Soft Convex Moderma Flex Pouch	12
2281055	One-Piece Closed – Soft Convex Moderma Flex Pouch	12
2281125	One-Piece Closed – Soft Convex Moderma Flex Pouch	12
2281130	One-Piece Closed – Soft Convex Moderma Flex Pouch	12
2281135	One-Piece Closed – Soft Convex Moderma Flex Pouch	12
2281140	One-Piece Closed – Soft Convex Moderma Flex Pouch	12
2909064	One-Piece Urostomy – Flat Moderma Flex Multi-Chamber Pouch – with Tape Border	28
2909125	One-Piece Urostomy – Flat Moderma Flex Multi-Chamber Pouch – with Tape Border	28
2909130	One-Piece Urostomy – Flat Moderma Flex Multi-Chamber Pouch – with Tape Border	28
2909135	One-Piece Urostomy – Flat Moderma Flex Multi-Chamber Pouch – with Tape Border	28
2949064	One-Piece Urostomy – Flat Moderma Flex Multi-Chamber Pouch – with Tape Border	28
2969025	One-Piece Urostomy – Firm Convex Moderma Flex Multi-Chamber Pouch – with Tape Border	32
2969038	One-Piece Urostomy – Firm Convex Moderma Flex Multi-Chamber Pouch – with Tape Border	32
2969120	One-Piece Urostomy – Firm Convex Moderma Flex Multi-Chamber Pouch – with Tape Border	32
2969125	One-Piece Urostomy – Firm Convex Moderma Flex Multi-Chamber Pouch – with Tape Border	32
2969130	One-Piece Urostomy – Firm Convex Moderma Flex Multi-Chamber Pouch – with Tape Border	32
2969135	One-Piece Urostomy – Firm Convex Moderma Flex Multi-Chamber Pouch – with Tape Border	32
2969140	One-Piece Urostomy – Firm Convex Moderma Flex Multi-Chamber Pouch – with Tape Border	32
2979025	One-Piece Urostomy – Firm Convex Moderma Flex Multi-Chamber Pouch – with Tape Border	32
2979038	One-Piece Urostomy – Firm Convex Moderma Flex Multi-Chamber Pouch – with Tape Border	32
2979120	One-Piece Urostomy – Firm Convex Moderma Flex Multi-Chamber Pouch – with Tape Border	32
2979125	One-Piece Urostomy – Firm Convex Moderma Flex Multi-Chamber Pouch – with Tape Border	32
2979130	One-Piece Urostomy – Firm Convex Moderma Flex Multi-Chamber Pouch – with Tape Border	32
2979135	One-Piece Urostomy – Firm Convex Moderma Flex Multi-Chamber Pouch – with Tape Border	32
2979140	One-Piece Urostomy – Firm Convex Moderma Flex Multi-Chamber Pouch – with Tape Border	32

Code	Name	Page No.
5681038	One-Piece Drainable – Soft Convex Moderma Flex Pouch	20
5681055	One-Piece Drainable – Soft Convex Moderma Flex Pouch	20
5681120	One-Piece Drainable – Soft Convex Moderma Flex Pouch	20
5681125	One-Piece Drainable – Soft Convex Moderma Flex Pouch	20
5681130	One-Piece Drainable – Soft Convex Moderma Flex Pouch	20
5681135	One-Piece Drainable – Soft Convex Moderma Flex Pouch	20

Dansac Code Index

Code	Name	Page No.
070-20	Accessories Dansac Seals	63
070-30	Accessories Dansac Seals	63
070-40	Accessories Dansac Seals	63
071-30	Accessories Dansac TRE – Seals	63
071-40	Accessories Dansac TRE – Seals	63
072-48	Accessories Dansac TRE – Seals	63
080-00	Accessories Dansac Nodor S	65
083-01	Accessories Dansac EasiSpray™ Adhesive Remover	65
090-02	Accessories Belts	65
090-03	Accessories Belts	65
725-20/30	Accessories Dansac GX-tra Seals	63
725-30/30	Accessories Dansac GX-tra Seals	63
725-40/30	Accessories Dansac GX-tra Seals	63
725-50/30	Accessories Dansac GX-tra Seals	63
801-20	One-Piece Closed – Flat Nova 1 Closed Flat Pouch	11
801-25	One-Piece Closed – Flat Nova 1 Closed Flat Pouch	11
801-30	One-Piece Closed – Flat Nova 1 Closed Flat Pouch	11
801-35	One-Piece Closed – Flat Nova 1 Closed Flat Pouch	11
801-40	One-Piece Closed – Flat Nova 1 Closed Flat Pouch	11
801-45	One-Piece Closed – Flat Nova 1 Closed Flat Pouch	11
801-50	One-Piece Closed – Flat Nova 1 Closed Flat Pouch	11
802-20	One-Piece Closed – Flat Nova 1 Closed Flat Pouch	11
802-25	One-Piece Closed – Flat Nova 1 Closed Flat Pouch	11
802-30	One-Piece Closed – Flat Nova 1 Closed Flat Pouch	11
802-35	One-Piece Closed – Flat Nova 1 Closed Flat Pouch	11
803-20	One-Piece Closed – Flat Nova 1 Closed Flat Pouch	11
803-25	One-Piece Closed – Flat Nova 1 Closed Flat Pouch	11
803-30	One-Piece Closed – Flat Nova 1 Closed Flat Pouch	11
803-35	One-Piece Closed – Flat Nova 1 Closed Flat Pouch	11
805-20	One-Piece Closed – Flat Nova 1 Closed Flat Pouch	11

Code	Name	Page No.
806-20	One-Piece Closed – Flat Nova 1 Closed Flat Pouch	11
815-15	One-Piece Drainable – Flat Nova 1 Foldup Flat Pouch	19
816-15	One-Piece Drainable – Flat Nova 1 Foldup Flat Pouch	19
818-00	Paediatric Pouches – SoftFlex Skin Barrier Nova 1 Infant/Paediatric	57
818-10	Paediatric Pouches – SoftFlex Skin Barrier Nova 1 Infant/Paediatric	57
823-15	One-Piece Drainable – Flat Nova 1 Foldup Flat Pouch	19
823-25	One-Piece Drainable – Flat Nova 1 Foldup Flat Pouch	19
823-30	One-Piece Drainable – Flat Nova 1 Foldup Flat Pouch	19
823-35	One-Piece Drainable – Flat Nova 1 Foldup Flat Pouch	19
823-40	One-Piece Drainable – Flat Nova 1 Foldup Flat Pouch	19
824-15	One-Piece Drainable – Flat Nova 1 Foldup Flat Pouch	19
825-15	One-Piece Drainable – Flat Nova 1 Foldup Flat Pouch	19
826-15	One-Piece Drainable – Flat Nova 1 Foldup Flat Pouch	19
829-20	One-Piece Closed – Flat Nova 1 Closed Flat Pouch	11
829-25	One-Piece Closed – Flat Nova 1 Closed Flat Pouch	11
829-30	One-Piece Closed – Flat Nova 1 Closed Flat Pouch	11
829-40	One-Piece Closed – Flat Nova 1 Closed Flat Pouch	11
831-25	One-Piece Closed – Firm Convex Nova 1 Closed Convex Pouch	15
831-30	One-Piece Closed – Firm Convex Nova 1 Closed Convex Pouch	15
831-35	One-Piece Closed – Firm Convex Nova 1 Closed Convex Pouch	15
831-37	One-Piece Closed – Firm Convex Nova 1 Closed Convex Pouch	15
831-40	One-Piece Closed – Firm Convex Nova 1 Closed Convex Pouch	15
831-45	One-Piece Closed – Firm Convex Nova 1 Closed Convex Pouch	15
831-46	One-Piece Closed – Firm Convex Nova 1 Closed Convex Pouch	15
832-37	One-Piece Closed – Firm Convex Nova 1 Closed Convex Pouch	15
841-20	One-Piece Drainable – Firm Convex Nova 1 FoldUp Convex Pouch	27
841-25	One-Piece Drainable – Firm Convex Nova 1 FoldUp Convex Pouch	27
841-30	One-Piece Drainable – Firm Convex Nova 1 FoldUp Convex Pouch	27
841-35	One-Piece Drainable – Firm Convex Nova 1 FoldUp Convex Pouch	27
841-37	One-Piece Drainable – Firm Convex Nova 1 FoldUp Convex Pouch	27
841-46	One-Piece Drainable – Firm Convex Nova 1 FoldUp Convex Pouch	27
842-25	One-Piece Drainable – Firm Convex Nova 1 FoldUp Convex Pouch	27
842-30	One-Piece Drainable – Firm Convex Nova 1 FoldUp Convex Pouch	27
842-35	One-Piece Drainable – Firm Convex Nova 1 FoldUp Convex Pouch	27
842-37	One-Piece Drainable – Firm Convex Nova 1 FoldUp Convex Pouch	27
842-46	One-Piece Drainable – Firm Convex Nova 1 FoldUp Convex Pouch	27
861-15	One-Piece Drainable – Flat Nova 1 Foldup X3 Pouch	19
861-30	One-Piece Drainable – Flat Nova 1 Foldup X3 Pouch	19
871-25	One-Piece Closed – Soft Convex Nova 1 Closed Soft Mouldable Pouch	13
871-35	One-Piece Closed – Soft Convex Nova 1 Closed Soft Mouldable Pouch	13
881-25	One-Piece Drainable – Soft Mouldable Nova 1 FoldUp Soft Mouldable Pouch	23
881-30	One-Piece Drainable – Soft Mouldable Nova 1 FoldUp Soft Mouldable Pouch	23

Code	Name	Page No.
881-35	One-Piece Drainable – Soft Mouldable Nova 1 FoldUp Soft Mouldable Pouch	23
891-12	One-Piece Urostomy Nova 1 Urostomy X3 Pouch	31
891-20	One-Piece Urostomy Nova 1 Urostomy X3 Pouch	31
891-25	One-Piece Urostomy Nova 1 Urostomy X3 Pouch	31
891-30	One-Piece Urostomy Nova 1 Urostomy X3 Pouch	31
894-25	One-Piece Urostomy – Firm Convex Nova 1 Urostomy Convex Pouch	33
894-30	One-Piece Urostomy – Firm Convex Nova 1 Urostomy Convex Pouch	33
894-37	One-Piece Urostomy – Firm Convex Nova 1 Urostomy Convex Pouch	33
894-46	One-Piece Urostomy – Firm Convex Nova 1 Urostomy Convex Pouch	33
902-10	One-Piece Closed – Flat NovaLife GX 1 Closed Flat Pouch	9
902-25	One-Piece Closed – Flat NovaLife GX 1 Closed Flat Pouch	9
902-30	One-Piece Closed – Flat NovaLife GX 1 Closed Flat Pouch	9
902-35	One-Piece Closed – Flat NovaLife GX 1 Closed Flat Pouch	9
902-40	One-Piece Closed – Flat NovaLife GX 1 Closed Flat Pouch	9
902-45	One-Piece Closed – Flat NovaLife GX 1 Closed Flat Pouch	9
902-50	One-Piece Closed – Flat NovaLife GX 1 Closed Flat Pouch	9
903-30	One-Piece Closed – Flat NovaLife GX 1 Closed Flat Pouch	9
903-35	One-Piece Closed – Flat NovaLife GX 1 Closed Flat Pouch	9
906-10	One-Piece Closed – Flat NovaLife GX 1 Closed Flat Pouch	9
915-10	One-Piece Drainable – Flat NovaLife GX 1 Open Flat Pouch	17
916-10	One-Piece Drainable – Flat NovaLife GX 1 Open Flat Pouch	17
923-10	One-Piece Drainable – Flat NovaLife GX 1 Open Flat Pouch	17
923-25	One-Piece Drainable – Flat NovaLife GX 1 Open Flat Pouch	17
923-30	One-Piece Drainable – Flat NovaLife GX 1 Open Flat Pouch	17
923-35	One-Piece Drainable – Flat NovaLife GX 1 Open Flat Pouch	17
924-10	One-Piece Drainable – Flat NovaLife GX 1 Open Flat Pouch	17
925-10	One-Piece Drainable – Flat NovaLife GX 1 Open Flat Pouch	17
932-25	One-Piece Closed – Firm Convex NovaLife GX 1 Closed Convex Pouch	15
932-30	One-Piece Closed – Firm Convex NovaLife GX 1 Closed Convex Pouch	15
932-35	One-Piece Closed – Firm Convex NovaLife GX 1 Closed Convex Pouch	15
932-40	One-Piece Closed – Firm Convex NovaLife GX 1 Closed Convex Pouch	15
932-45	One-Piece Closed – Firm Convex NovaLife GX 1 Closed Convex Pouch	15
941-25	One-Piece Drainable – Firm Convex NovaLife GX 1 Open Convex Pouch	25
941-30	One-Piece Drainable – Firm Convex NovaLife GX 1 Open Convex Pouch	25
941-35	One-Piece Drainable – Firm Convex NovaLife GX 1 Open Convex Pouch	25
941-37	One-Piece Drainable – Firm Convex NovaLife GX 1 Open Convex Pouch	25
941-46	One-Piece Drainable – Firm Convex NovaLife GX 1 Open Convex Pouch	25
941-50	One-Piece Drainable – Firm Convex NovaLife GX 1 Open Convex Pouch	25
941-59	One-Piece Drainable – Firm Convex NovaLife GX 1 Open Convex Pouch	25
942-30	One-Piece Drainable – Firm Convex NovaLife GX 1 Open Convex Pouch	25
942-35	One-Piece Drainable – Firm Convex NovaLife GX 1 Open Convex Pouch	25
942-37	One-Piece Drainable – Firm Convex NovaLife GX 1 Open Convex Pouch	25

Code	Name	Page No.
942-46	One-Piece Drainable – Firm Convex NovaLife GX 1 Open Convex Pouch	25
942-59	One-Piece Drainable – Firm Convex NovaLife GX 1 Open Convex Pouch	25
971-25	One-Piece Closed – Soft Convex NovaLife GX 1 Closed Soft Convex Pouch	13
971-30	One-Piece Closed – Soft Convex NovaLife GX 1 Closed Soft Convex Pouch	13
971-34	One-Piece Closed – Soft Convex NovaLife GX 1 Closed Soft Convex Pouch	13
971-35	One-Piece Closed – Soft Convex NovaLife GX 1 Closed Soft Convex Pouch	13
971-40	One-Piece Closed – Soft Convex NovaLife GX 1 Closed Soft Convex Pouch	13
971-44	One-Piece Closed – Soft Convex NovaLife GX 1 Closed Soft Convex Pouch	13
971-54	One-Piece Closed – Soft Convex NovaLife GX 1 Closed Soft Convex Pouch	13
981-24	One-Piece Drainable – Soft Convex NovaLife GX 1 Open Soft Convex Pouch	21
981-25	One-Piece Drainable – Soft Convex NovaLife GX 1 Open Soft Convex Pouch	21
981-30	One-Piece Drainable – Soft Convex NovaLife GX 1 Open Soft Convex Pouch	21
981-34	One-Piece Drainable – Soft Convex NovaLife GX 1 Open Soft Convex Pouch	21
981-35	One-Piece Drainable – Soft Convex NovaLife GX 1 Open Soft Convex Pouch	21
981-40	One-Piece Drainable – Soft Convex NovaLife GX 1 Open Soft Convex Pouch	21
981-44	One-Piece Drainable – Soft Convex NovaLife GX 1 Open Soft Convex Pouch	21
981-45	One-Piece Drainable – Soft Convex NovaLife GX 1 Open Soft Convex Pouch	21
981-54	One-Piece Drainable – Soft Convex NovaLife GX 1 Open Soft Convex Pouch	21
983-24	One-Piece Drainable – Soft Convex NovaLife GX 1 Open Soft Convex Pouch	21
983-30	One-Piece Drainable – Soft Convex NovaLife GX 1 Open Soft Convex Pouch	21
983-34	One-Piece Drainable – Soft Convex NovaLife GX 1 Open Soft Convex Pouch	21
983-35	One-Piece Drainable – Soft Convex NovaLife GX 1 Open Soft Convex Pouch	21
1136-20	Two-Piece Barriers – Flat Nova 2 Flat Wafer	39
1143-15	Two-Piece Barriers – Flat Nova 2 Flat Wafer	39
1143-25	Two-Piece Barriers – Flat Nova 2 Flat Wafer	39
1143-30	Two-Piece Barriers – Flat Nova 2 Flat Wafer	39
1155-15	Two-Piece Barriers – Flat Nova 2 Flat Wafer	39
1155-35	Two-Piece Barriers – Flat Nova 2 Flat Wafer	39
1155-40	Two-Piece Barriers – Flat Nova 2 Flat Wafer	39
1170-15	Two-Piece Barriers – Flat Nova 2 Flat Wafer	39
1170-45	Two-Piece Barriers – Flat Nova 2 Flat Wafer	39
1170-50	Two-Piece Barriers – Flat Nova 2 Flat Wafer	39
1201-36	Two-Piece Pouches – Closed Nova 2 Closed Pouch	49
1201-43	Two-Piece Pouches – Closed Nova 2 Closed Pouch	49
1201-55	Two-Piece Pouches – Closed Nova 2 Closed Pouch	49
1201-70	Two-Piece Pouches – Closed Nova 2 Closed Pouch	49
1203-36	Two-Piece Pouches – Closed Nova 2 Closed Pouch	49
1215-43	Two-Piece Pouches – Drainable Nova 2 FoldUp Pouch	51
1215-55	Two-Piece Pouches – Drainable Nova 2 FoldUp Pouch	51
1215-70	Two-Piece Pouches – Drainable Nova 2 FoldUp Pouch	51
1217-43	Two-Piece Pouches – Urostomy Nova 2 Urostomy Pouch	53
1217-55	Two-Piece Pouches – Urostomy Nova 2 Urostomy Pouch	53

Code	Name	Page No.
1218-43	Two-Piece Pouches – Urostomy Nova 2 Urostomy Pouch	53
1218-55	Two-Piece Pouches – Urostomy Nova 2 Urostomy Pouch	53
1301-36	Two-Piece Pouches – Closed NovaLife 2 Closed Pouch	49
1301-43	Two-Piece Pouches – Closed NovaLife 2 Closed Pouch	49
1301-55	Two-Piece Pouches – Closed NovaLife 2 Closed Pouch	49
1301-70	Two-Piece Pouches – Closed NovaLife 2 Closed Pouch	49
1303-43	Two-Piece Pouches – Closed NovaLife 2 Closed Pouch	49
1315-43	Two-Piece Pouches – Drainable NovaLife 2 Open Pouch	51
1315-55	Two-Piece Pouches – Drainable NovaLife 2 Open Pouch	51
1315-70	Two-Piece Pouches – Drainable NovaLife 2 Open Pouch	51
1321-43	Two-Piece Pouches – Drainable NovaLife 2 Open Pouch	51
1321-55	Two-Piece Pouches – Drainable NovaLife 2 Open Pouch	51
1343-15	Two-Piece Barriers Nova 2 X3 Wafer	41
1343-25	Two-Piece Barriers Nova 2 X3 Wafer	41
1355-15	Two-Piece Barriers Nova 2 X3 Wafer	41
1355-35	Two-Piece Barriers Nova 2 X3 Wafer	41
1443-21	Two-Piece Barriers Nova 2 Soft Mouldable Wafer	41
1443-25	Two-Piece Barriers Nova 2 Soft Mouldable Wafer	41
1443-28	Two-Piece Barriers Nova 2 Soft Mouldable Wafer	41
1455-32	Two-Piece Barriers Nova 2 Soft Mouldable Wafer	41
1455-38	Two-Piece Barriers Nova 2 Soft Mouldable Wafer	41
1543-15	Two-Piece Barriers – Firm Convex Nova 2 Convex Wafer	43
1543-25	Two-Piece Barriers – Firm Convex Nova 2 Convex Wafer	43
1543-30	Two-Piece Barriers – Firm Convex Nova 2 Convex Wafer	43
1555-15	Two-Piece Barriers – Firm Convex Nova 2 Convex Wafer	43
1555-35	Two-Piece Barriers – Firm Convex Nova 2 Convex Wafer	43
1555-40	Two-Piece Barriers – Firm Convex Nova 2 Convex Wafer	43
1743-25	Two-Piece Barriers – Flat NovaLife GX 2 Flat Wafer	39
1755-35	Two-Piece Barriers – Flat NovaLife GX 2 Flat Wafer	39
1770-10	Two-Piece Barriers – Flat NovaLife GX 2 Flat Wafer	39
1770-50	Two-Piece Barriers – Flat NovaLife GX 2 Flat Wafer	39
1936-20	Two-Piece Barriers – Firm Convex NovaLife GX 2 Convex Wafer	43
1936-23	Two-Piece Barriers – Firm Convex NovaLife GX 2 Convex Wafer	43
1943-10	Two-Piece Barriers – Firm Convex NovaLife GX 2 Convex Wafer	43
1943-25	Two-Piece Barriers – Firm Convex NovaLife GX 2 Convex Wafer	43
1943-30	Two-Piece Barriers – Firm Convex NovaLife GX 2 Convex Wafer	43
1955-10	Two-Piece Barriers – Firm Convex NovaLife GX 2 Convex Wafer	43
1955-35	Two-Piece Barriers – Firm Convex NovaLife GX 2 Convex Wafer	43
1955-40	Two-Piece Barriers – Firm Convex NovaLife GX 2 Convex Wafer	43
1970-10	Two-Piece Barriers – Firm Convex NovaLife GX 2 Convex Wafer	43
1970-45	Two-Piece Barriers – Firm Convex NovaLife GX 2 Convex Wafer	43
3001-15	One-Piece Closed – Flat NovaLife TRE 1 Closed Flat Pouch	9

Code	Name	Page No.
3001-25	One-Piece Closed – Flat NovaLife TRE 1 Closed Flat Pouch	9
3001-30	One-Piece Closed – Flat NovaLife TRE 1 Closed Flat Pouch	9
3001-35	One-Piece Closed – Flat NovaLife TRE 1 Closed Flat Pouch	9
3001-40	One-Piece Closed – Flat NovaLife TRE 1 Closed Flat Pouch	9
3001-45	One-Piece Closed – Flat NovaLife TRE 1 Closed Flat Pouch	9
3005-15	One-Piece Closed – Flat NovaLife TRE 1 Closed Flat Pouch	9
3005-25	One-Piece Closed – Flat NovaLife TRE 1 Closed Flat Pouch	9
3005-30	One-Piece Closed – Flat NovaLife TRE 1 Closed Flat Pouch	9
3007-15	One-Piece Closed – Flat NovaLife TRE 1 Closed Flat Pouch	9
3007-25	One-Piece Closed – Flat NovaLife TRE 1 Closed Flat Pouch	9
3007-30	One-Piece Closed – Flat NovaLife TRE 1 Closed Flat Pouch	9
3015-15	One-Piece Drainable – Flat NovaLife TRE 1 Open Flat Pouch	17
3015-25	One-Piece Drainable – Flat NovaLife TRE 1 Open Flat Pouch	17
3015-30	One-Piece Drainable – Flat NovaLife TRE 1 Open Flat Pouch	17
3016-15	One-Piece Drainable – Flat NovaLife TRE 1 Open Flat Pouch	17
3017-15	One-Piece Drainable – Flat NovaLife TRE 1 Open Flat Pouch	17
3017-25	One-Piece Drainable – Flat NovaLife TRE 1 Open Flat Pouch	17
3017-30	One-Piece Drainable – Flat NovaLife TRE 1 Open Flat Pouch	17
3023-15	One-Piece Drainable – Flat NovaLife TRE 1 Open Flat Pouch	17
3023-25	One-Piece Drainable – Flat NovaLife TRE 1 Open Flat Pouch	17
3023-30	One-Piece Drainable – Flat NovaLife TRE 1 Open Flat Pouch	17
3023-35	One-Piece Drainable – Flat NovaLife TRE 1 Open Flat Pouch	17
3023-40	One-Piece Drainable – Flat NovaLife TRE 1 Open Flat Pouch	17
3024-15	One-Piece Drainable – Flat NovaLife TRE 1 Open Flat Pouch	17
3025-15	One-Piece Drainable – Flat NovaLife TRE 1 Open Flat Pouch	17
3026-15	One-Piece Drainable – Flat NovaLife TRE 1 Open Flat Pouch	17
3031-30	One-Piece Closed – Firm Convex NovaLife TRE 1 Closed Convex Pouch	15
3031-35	One-Piece Closed – Firm Convex NovaLife TRE 1 Closed Convex Pouch	15
3031-37	One-Piece Closed – Firm Convex NovaLife TRE 1 Closed Convex Pouch	15
3031-59	One-Piece Closed – Firm Convex NovaLife TRE 1 Closed Convex Pouch	15
3035-25	One-Piece Closed – Firm Convex NovaLife TRE 1 Closed Convex Pouch	15
3041-30	One-Piece Drainable – Firm Convex NovaLife TRE 1 Open Convex Pouch	25
3041-35	One-Piece Drainable – Firm Convex NovaLife TRE 1 Open Convex Pouch	25
3041-37	One-Piece Drainable – Firm Convex NovaLife TRE 1 Open Convex Pouch	25
3041-40	One-Piece Drainable – Firm Convex NovaLife TRE 1 Open Convex Pouch	25
3041-46	One-Piece Drainable – Firm Convex NovaLife TRE 1 Open Convex Pouch	25
3041-59	One-Piece Drainable – Firm Convex NovaLife TRE 1 Open Convex Pouch	25
3042-37	One-Piece Drainable – Firm Convex NovaLife TRE 1 Open Convex Pouch	25
3042-59	One-Piece Drainable – Firm Convex NovaLife TRE 1 Open Convex Pouch	25
3045-25	One-Piece Drainable – Firm Convex NovaLife TRE 1 Open Convex Pouch	25
3071-24	One-Piece Closed – Soft Convex NovaLife TRE 1 Closed Soft Convex Pouch	13
3071-25	One-Piece Closed – Soft Convex NovaLife TRE 1 Closed Soft Convex Pouch	13

Code	Name	Page No.
3071-30	One-Piece Closed – Soft Convex NovaLife TRE 1 Closed Soft Convex Pouch	13
3071-34	One-Piece Closed – Soft Convex NovaLife TRE 1 Closed Soft Convex Pouch	13
3071-35	One-Piece Closed – Soft Convex NovaLife TRE 1 Closed Soft Convex Pouch	13
3071-40	One-Piece Closed – Soft Convex NovaLife TRE 1 Closed Soft Convex Pouch	13
3071-44	One-Piece Closed – Soft Convex NovaLife TRE 1 Closed Soft Convex Pouch	13
3071-54	One-Piece Closed – Soft Convex NovaLife TRE 1 Closed Soft Convex Pouch	13
3081-24	One-Piece Drainable – Soft Convex NovaLife TRE 1 Open Soft Convex Pouch	21
3081-25	One-Piece Drainable – Soft Convex NovaLife TRE 1 Open Soft Convex Pouch	21
3081-30	One-Piece Drainable – Soft Convex NovaLife TRE 1 Open Soft Convex Pouch	21
3081-34	One-Piece Drainable – Soft Convex NovaLife TRE 1 Open Soft Convex Pouch	21
3081-35	One-Piece Drainable – Soft Convex NovaLife TRE 1 Open Soft Convex Pouch	21
3081-40	One-Piece Drainable – Soft Convex NovaLife TRE 1 Open Soft Convex Pouch	21
3081-44	One-Piece Drainable – Soft Convex NovaLife TRE 1 Open Soft Convex Pouch	21
3081-45	One-Piece Drainable – Soft Convex NovaLife TRE 1 Open Soft Convex Pouch	21
3081-54	One-Piece Drainable – Soft Convex NovaLife TRE 1 Open Soft Convex Pouch	21
3083-25	One-Piece Drainable – Soft Convex NovaLife TRE 1 Open Soft Convex Pouch	21
3083-30	One-Piece Drainable – Soft Convex NovaLife TRE 1 Open Soft Convex Pouch	21
3083-34	One-Piece Drainable – Soft Convex NovaLife TRE 1 Open Soft Convex Pouch	21
3083-35	One-Piece Drainable – Soft Convex NovaLife TRE 1 Open Soft Convex Pouch	21
3084-34	One-Piece Drainable – Soft Convex NovaLife TRE 1 Open Soft Convex Pouch	21
3143-15	Two-Piece Barriers – Flat NovaLife TRE 2 Flat Wafer	39
3143-25	Two-Piece Barriers – Flat NovaLife TRE 2 Flat Wafer	39
3143-30	Two-Piece Barriers – Flat NovaLife TRE 2 Flat Wafer	39
3155-15	Two-Piece Barriers – Flat NovaLife TRE 2 Flat Wafer	39
3243-15	Two-Piece Barriers – Firm Convex NovaLife TRE 2 Convex Wafer	43
3243-25	Two-Piece Barriers – Firm Convex NovaLife TRE 2 Convex Wafer	43
3243-30	Two-Piece Barriers – Firm Convex NovaLife TRE 2 Convex Wafer	43
3255-15	Two-Piece Barriers – Firm Convex NovaLife TRE 2 Convex Wafer	43
3255-35	Two-Piece Barriers – Firm Convex NovaLife TRE 2 Convex Wafer	43
3255-40	Two-Piece Barriers – Firm Convex NovaLife TRE 2 Convex Wafer	43
9010-10	One-Piece Closed – Flat NovaLife GX 1 Closed Flat Pouch	9
9010-25	One-Piece Closed – Flat NovaLife GX 1 Closed Flat Pouch	9
9010-30	One-Piece Closed – Flat NovaLife GX 1 Closed Flat Pouch	9
9010-35	One-Piece Closed – Flat NovaLife GX 1 Closed Flat Pouch	9
9010-40	One-Piece Closed – Flat NovaLife GX 1 Closed Flat Pouch	9
9010-45	One-Piece Closed – Flat NovaLife GX 1 Closed Flat Pouch	9
9010-50	One-Piece Closed – Flat NovaLife GX 1 Closed Flat Pouch	9
9050-10	One-Piece Closed – Flat NovaLife GX 1 Closed Flat Pouch	9
9050-25	One-Piece Closed – Flat NovaLife GX 1 Closed Flat Pouch	9
9050-30	One-Piece Closed – Flat NovaLife GX 1 Closed Flat Pouch	9
9050-35	One-Piece Closed – Flat NovaLife GX 1 Closed Flat Pouch	9
9310-25	One-Piece Closed – Firm Convex NovaLife GX 1 Closed Convex Pouch	15

Code	Name	Page No.
9310-30	One-Piece Closed – Firm Convex NovaLife GX 1 Closed Convex Pouch	15
9310-35	One-Piece Closed – Firm Convex NovaLife GX 1 Closed Convex Pouch	15
9310-37	One-Piece Closed – Firm Convex NovaLife GX 1 Closed Convex Pouch	15
9310-45	One-Piece Closed – Firm Convex NovaLife GX 1 Closed Convex Pouch	15
9310-50	One-Piece Closed – Firm Convex NovaLife GX 1 Closed Convex Pouch	15
9510-25	One-Piece Drainable – Firm Convex NovaLife GX+ 1 Open Convex Pouch	27
9510-30	One-Piece Drainable – Firm Convex NovaLife GX+ 1 Open Convex Pouch	27
9510-37	One-Piece Drainable – Firm Convex NovaLife GX+ 1 Open Convex Pouch	27
9510-59	One-Piece Drainable – Firm Convex NovaLife GX+ 1 Open Convex Pouch	27
9610-10	One-Piece Drainable – Flat NovaLife GX+ 1 Open Flat Pouch	19
9610-25	One-Piece Drainable – Flat NovaLife GX+ 1 Open Flat Pouch	19
9630-25	One-Piece Drainable – Flat NovaLife GX+ 1 Open Flat Pouch	19
9630-30	One-Piece Drainable – Flat NovaLife GX+ 1 Open Flat Pouch	19
9630-35	One-Piece Drainable – Flat NovaLife GX+ 1 Open Flat Pouch	19
70000-0000	Accessories Dansac Skin Lotion	63
71000-0000	Accessories Dansac Skin Lotion Wipes	65
77550-0	Accessories Dansac Soft Paste	63
W9759	Accessories Windless Gas Suppressant	65
W9760	Accessories Windless Gas Suppressant	65

Prior to use, be sure to read the Instructions for Use for information regarding Intended Use, Contraindications, Warnings, Precautions and Instructions.

The Dansac logo, NovaLife, TRE, and Stoma, Skin, Soul are trademarks of Dansac A/S. All other trademarks and copyrights are the property of their respective owners. © 2020 Dansac A/S.

The Hollister Logo, Adapt, Apogee, CeraPlus and "Healthy skin. Positive outcomes." are trademarks of Hollister Incorporated. CenterPointLock coupling system (trademark) for all base plates. All other trademarks and copyrights are the property of their respective owners. Not all products are CE marked. © 2020 Hollister Incorporated.

Printed in Australia. 2020/06.